

MRI Abdomen/Pelvis W & WO												
Technical Parameters												
Anatomic Coverage:	Entire liver through perineum											
Coil:	Torso											
Scanner Preference:	1.5T											
Contrast:	Pre-injection:								Post-injection: trigger when bolus reached SMA			
Scanning Plane:	Cor Abd/Pel	Ax Abd/Pel	Ax Abd	Ax Pel	Ax Abd/Pel	Ax Abd/Pel	Ax Pel	Ax Abd	Ax Abd	Ax Pel	Ax Abd (4m delay)	Sag Pel
Slice Thickness:	7mm	7mm	7mm	7mm	5mm	7mm	7mm	7mm	7mm	7mm	7mm	5mm
Gap:	1mm	1mm	1mm	1mm	1mm	1mm	1mm	1mm	1mm	1mm	1mm	1mm
Weighting:	T2	T1	T2	T2	T2	T2	T1	T1	T1	T1	T1	T2
Mode:	SS BH	Dual Echo SPGR BH	DWI	DWI	SS	TSE RT	3D THRIVE BH	3D mDixon PRE BH	3D mDixon x3 phases BH	3D THRIVE BH	3D mDixon BH	TSE
FAT SAT:	NO	NO	SPIR	SPIR	NO	SPAIR	YES	YES	YES	YES	YES	NO
FOV:	Fit to patient	Fit to patient	Fit to patient	Fit to patient	20-24mm	Fit to patient	Fit to patient	Fit to patient	Fit to patient	Fit to patient	Fit to patient	24mm
MPR:	NO	NO	NO	NO	NO	NO	NO	NO	NO	NO	NO	NO
Notes:												