

MRI Female Pelvis W & WO												
Technical Parameters												
Anatomic Coverage:	Iliac crest to perinum											
Coil:	Torso											
Scanner Preference:	1.5T											
Contrast:	Pre-injection:								Post-injection:			
Scanning Plane:	Coronal	Axial	Axial	Sagittal	Coronal	Axial PRE	Axial	Axial	Axial POST	Sagittal	Coronal	Axial
Slice Thickness:	5mm	5mm	5mm	5mm	5mm	5mm	5mm	5mm	5mm	5mm	5mm	5mm
Gap:	1mm	1mm	1mm	1mm	1mm	1mm	1mm	1mm	1mm	1mm	1mm	1mm
Weighting:	T2	T1	T2	T2	T2	T1	T2	T2	T1	T1	T1	T2
Mode:	SS-TSE BH	TSE	TSE	TSE	TSE FS	3D THRIVE	DWI/ADC	IN/OUT Phase	3D THRIVE High Res	3D THRIVE High Res	3D THRIVE High Res	DWI
FAT SAT:	No	NO	No	NO	YES	YES	NO	NO	YES	YES	YES	SPIR
FOV:	Top of kidneys to pelvis	20-24cm	20-24cm	24cm	20-24cm	Match AX TSE T2 Pelvis	Match AX TSE T2 Pelvis	Match AX TSE T2 Pelvis	Match high res precontrast THRIVE	Match SAG TSE T2 Pelvis	Match COR TSE FS T2 Pelvis	Match AX TSE T2 Pelvis
MPR:												
Notes:							*If indication is for endometriosis or ovairan mass					