

MRI Pancreas & MRCP W & WO												
Technical Parameters												
Anatomic Coverage:	Diaphragm to aortic bifurcation											
Coil:	Torso											
Scanner Preference:	1.5T or 3T											
Contrast:	Pre-injection:									Post-injection: trigger when bolus reaches SMA		
Scanning Plane:	Coronal	Axial	Axial	Axial	Axial	Coronal	Coronal	Axial	Axial	Axial	Coronal	Axial
Slice Thickness:	7mm	7mm	4mm	7mm	5mm	1.4mm	12mm	7mm	5mm	5mm	5mm	5mm
Gap:	1mm	1mm	-2mm	1mm	0.5mm	0mm		1mm	0.5mm	0.5mm	0.5mm	0.5mm
Weighting:	T2	T1	T2	T1	T2	T2	T2	T2	T2	T1	T1	T1
Mode:	SSFSE BH	Dual Echo SPGR BH	bTFE BH	3D mDixon pre BH	TSE RT	3D MRCP	Radial MRCP	DWI	SSFSE RT	3D mDixon x3 phases BH	3D mDixon 3min post	3D mDixon BH
FAT SAT:	NO	NO	NO	YES	SPAIR	YES	YES	SPIR	YES	YES	YES	YES
FOV:	Liver through bottom of kidneys	Entire Liver		Entire Liver	Entire Liver	Bottom 2/3 of liver through bottom of pancreas	Liver through bottom of kidneys	Liver through bottom of kidneys	Entire Liver	Entire Liver	Diaphragm to aortic bifurcation	Entire Liver
MPR:	NO	NO	NO	NO	NO	NO	NO	NO	NO	NO	NO	NO
Notes:												