



Ultra-Light

EXCELLENT PERFORMANCE & MANEUVERABILITY

The Küschall Ultra-Light is a traditional cross folding wheelchair, offering you a more active appearance with its fixed front frame. Since it's first release 20 years ago, the Ultra-Light has evolved to keep up with the pace of your busy modern day lifestyle.



Yes, you can.®

küschall®



Ultra-Light

One of the key characteristics of the Ultra Light is the effortless folding; this is possible due to the unique shape and design of the crossbrace. Specially positioned pivot points allows a smooth folding and unfolding action designed for people with limited strength; in addition, the built in rigidity makes it very light and energy efficient to propel.

The KÜSCHALL Ultra-Light has an optional foldable backrest which sits neatly over the seat rails which, when folded makes it easy to handle, move and transport. When folded, it fits into the smallest of storage spaces and car boots.

Features



Available with either 70 or 80° front frame angle.



Cross-foldable ultra-light weight wheelchair.



Only 29cm width when folded.



Rear wheel fixation.

Tech Specs



Seat width
34 - 44cm



Seat depth
32 - 50cm



Seat height
front: 41 - 53cm
rear: 38 - 50cm



Total width
seat width + 18cm



Total width folded
29cm



Total length
88cm



Max. user weight
130 kg



Transport weight
w/out rear wheels
7.5 kg



Total weight
unloaded
9.9kg



Turning radius
92cm



Max. safe slope
7°



Crash Test Approved
Fulfills the requirements of ISO 7176-19



FOR FURTHER INFORMATION:

New Zealand: Freephone 0800 468 222, email sales@invacare.co.nz or visit www.invacare.co.nz

Australia: Freephone 1800 460 460, email sales@invacare.com.au or visit www.invacare.com.au

Asia: Phone +61 2 8839 5330, email exportsales@invacare.com.au or visit www.invacare.com

Thailand: Phone +66 0 2108 8555 or visit www.invacare.com



Yes, you can.®