

ASO Materials Provided

MATERIALS	25 TEST KIT	50 TEST KIT	100 TEST KIT	1000 TEST KIT
ASO Latex Reagent	1.0 ml	2.0 ml	2 x 2.0 ml	20 x 2.0 ml
Reactive Control	0.5 ml	1.0 ml	1.0 ml	10 x 1.0 ml
Nonreactive Control	0.5 ml	1.0 ml	1.0 ml	10 x 1.0 ml
0.03 ml Disposable Stirrer Pipets	25	50	100	1000
Disposable Test Cards (6-well)	5	9	17	170

ADDITIONAL MATERIALS REQUIRED

- + 13 x 75 mm test tubes
- + Timer
- + Volumetric pipet to deliver 0.25 ml
- + Saline (0.9% NaCl solution)
- + Mechanical rotator (optional)