

RF Materials Provided

MATERIALS	25 TEST KIT	100 TEST KIT	1000 TEST KIT
RF Direct Latex Reagent	1.0 ml	4.0 ml	10 x 4.0 ml
Reactive Control	0.5 ml	1.0 ml	10 x 1.0 ml
Nonreactive Control	0.5 ml	1.0 ml	10 x 1.0 ml
Glycine-Saline Solution	50 ml	100 ml	10 x 100 ml
0.05 ml Disposable Stirrer Pipets	25	100	1000
Disposable Test Cards (6-well)	5	17	170

ADDITIONAL MATERIALS REQUIRED

- + Volumetric pipet to deliver 0.25 ml
- + Timer
- + Mechanical rotator (optional)
- + 13 x 75 mm test tubes