

Product Catalog

Medical Solutions



ECG Electrodes



Defibrillation Electrodes



Neutral Electrodes
and ESU accessories



Content

About Leonhard Lang	04
■ ECG Electrodes	06
Stress Test & Holter	08
Resting ECG	09
Monitoring	10
Universal	11
MR-conditional & Radiotranslucent	13
Pediatrics & Neonates	14
■ Defibrillation Electrodes	16
Multifunctional Radiotranslucent Defibrillation Electrodes for Adults	18
Multifunctional Defibrillation Electrodes for Pediatrics	19
■ Neutral Electrodes	20
Neutral Electrodes	22
Neutral Electrodes with cables	22
Pencils & Tip Cleaner	23
Reusable Cables	23

About

Leonhard Lang's internationally renowned innovation and production center is located in the heart of the Tyrolean Alps.

Founded in 1872, Leonhard Lang has been operating as a family business for more than 150 years. Today, the focus is on medical technology. Over the years, joint top performance has resulted in production, special machine construction, process management, and research and development at the highest level. The electrode systems developed and produced by Leonhard Lang are in use worldwide.



We will always take the sustainable path


A resource-conserving approach to nature is important. The way we treat our environment today, determines our tomorrow. That is why Leonhard Lang takes action in all its business areas – for the health of future generations.

ECG Electrodes

Excellent quality ECG recordings are dependent on three main factors: good skin preparation, superior quality electrodes, and a continuous unbroken electrical pathway from the skin to the ECG machine.

Electrodes often have to perform under very strenuous conditions, therefore at Leonhard Lang we dedicate valuable time and attention to the design and construction of our electrodes; we then carefully select only the finest quality materials and component parts to produce them.

Uniquely, we design and build our own machines that are state of the art in manufacturing technology ensuring you quality products that are second to none. "Connect with quality" is our motto and our goal is to provide you with consistently reliable products.



Conductive Gels

Different ECG applications require gels that are precisely tailored to their specific performance demands. To ensure optimal performance for every application, Leonhard Lang develops and manufactures its conductive hydrogels in-house.

This guarantees maximum freshness, consistent quality and reliable conductivity. Our gels combine excellent signal response with skin-friendliness – supported by continuous innovation and full process control at our production site.

■ Aqua-Wet Gel = LIQUID GEL

Our Aqua-Wet Gel ensures particularly fast signal pick-up at the start of the ECG, regardless of skin preparation. Skintact® ECG electrodes with Aqua-Wet Gel are the best choice for short-term applications.

■ Aqua-Tac Gel = SOLID GEL

Skintact® ECG electrodes with Aqua-Tac Gel are used for both long-term and short-term applications. They can be removed from the skin without leaving any residue.

■ Aqua-Set Gel = SOLID WET GEL

Our Aqua-Set Gel combines the advantages of liquid and solid gel. Its special formula guarantees excellent electrical conductivity. The water content of Aqua-Set Gel ensures rapid reduction of skin resistance.

Backing Material

■ PE Foam

Demanding conditions, such as use in emergency care or for stress tests challenge the stability and reliability of ECG electrodes. Leonhard Lang's PE foams are the right choice – they are impermeable to liquids and thus protecting the sensor and gel. The thin PE material allows maximum flexibility.

■ Soft Cloth Backing

Skintact® ECG electrodes made of soft cloth backing offer all the advantages of nonwoven electrodes, but with significantly increased elasticity. Thus, they adapt well to body movements. Electrodes made of soft cloth backing offer a high level of wearing comfort for long-term applications and adhere reliably even to active patients.

■ Microporous Tape

The skin-friendly, microporous tape adapts well to the contours of the body and thus increases wearing comfort. It also allows moisture to pass through. Skintact® ECG electrodes with microporous tape allow the skin to breathe and transpire – which is particularly important for long-term applications.

■ Clear Tape

The clear tape is transparent. The patient's skin can therefore be observed during the entire wearing period. Clear tape reliably protects the sensor area against liquids.

Connectors

■ Stud

■ Tab

■ Tab'n'Snap

■ 4 mm Banana Adapter

Information

All our Skintact® products are tested according to applicable standards, ensuring maximum safety and reliable performance.

Skintact® ECG electrodes comply with standard ANSI/AAMI EC12 for ECG electrodes, ISO 10993 for biocompatibility, and the European Medical Device Regulation (MDR).

Stress Test & Holter

■ Bonding safe

Skintact® stress test electrodes with reinforced adhesive stick reliably. Even when patients transpire heavily. The PE foam backing is also impermeable to liquids. The ECG electrodes are therefore also suitable for echocardiograms.

■ Comfortable to wear

The particularly skin-friendly Skintact® ECG electrodes T-604 B, T-VO01 and T-VO02 allow the skin to breathe and provide excellent wearing comfort, even during prolonged use.

■ Easy to handle


The easy-to-open pouch and the peel-off aid of each Skintact® ECG electrode facilitate daily handling. On the Skintact® offset electrodes, the cables can be connected even more conveniently and reduce the risk of squeezing gel onto the electrode surface.

The PVC-free variant

Like all Skintact® ECG electrodes without pre-attached leads, also the Skintact® offset electrodes are PVC-free

* Innovative Technology


The adhesive ring of the T-VO01 and T-VO02, made of PE foam, protects the gel area and provides additional adhesive strength. In combination with the offset connection, it significantly reduces motion artifacts.



FS-VB01

Use	stress test
Backing	PE foam, enhanced adhesion
Dimensions	52 × 35 mm
Gel	liquid
Packaging	30/pouch, 1200/box, 4800/case


Premium electrode with offset connection



T-VO01 *

Use	holter, event recording
Backing	microporous tape
Dimensions	66 × 60 mm
Gel	liquid
Packaging	30/pouch, 300/box, 1800/case

Premium electrode with offset connection and with foam ring around gel surface; breathable material



FS-VB04

Use	stress test
Backing	PE foam, enhanced adhesion
Dimensions	52 × 35 mm
Gel	liquid
Packaging	30/pouch, 750/box, 3000/case


Premium electrode with offset connection and adapter for 4 mm banana plug



T-VO02 *

Use	holter, event recording
Backing	microporous tape
Dimensions	Ø 72 mm
Gel	liquid
Packaging	30/pouch, 300/box, 1800/case


Premium electrode with offset connection and with foam ring around gel surface; breathable material



FS-50

Use	stress test, emergency care, echocardiogram
Backing	PE foam, enhanced adhesion
Dimensions	Ø 50 mm
Gel	liquid
Packaging	30/pouch, 1200/box, 4800/case

Impermeable to liquids



T-604 B

Use	holter, long-term monitoring
Backing	microporous tape
Dimensions	Ø 60 mm
Gel	solid wet
Packaging	30/pouch, 1200/box, 4800/case

Breathable backing material, textured liner for skin preparation

Resting ECG

■ Excellent Conductivity

Skintact® Easitabs with Aqua-Tac gel provide reliable signal quality. The premium versions FS-WB00 and FS-WA00 feature the Aqua-Wet gel and are particularly suitable for resting ECG examinations, even under demanding conditions.

■ Bonding safe and convenient


All Skintact® Easitabs feel soft and stick absolutely reliably. Nevertheless, they are removed after application without irritation or loss of hair.

■ Superior Safety

Proper use of Skintact® Easitabs eliminates cross infection whereas suction electrodes can breed bacteria and fungi.

X-ray translucent and MR conditional


All Skintact® Easitabs can be safely used during imaging procedures.



RT14

Backing	Paper/PET
Dimensions	34 × 14 mm
Gel	solid
Packaging	100/pouch, 500/box, 5000/case

X-ray / MR conditional, Easitab, 10 pcs./card



FS-WB00 / 10

Backing	PE foam, enhanced adhesion
Dimensions	37 × 28 mm
Gel	liquid
Packaging	100/pouch, 500/box, 3000/case


Premium Easitab, particularly for demanding conditions, 10 pcs./card, X-ray / MR conditional



RT38

Backing	Paper/PET
Dimensions	34 × 20 mm
Gel	solid
Packaging	100/pouch, 500/box, 5000/case

Easitab, 10 pcs./card, X-ray / MR conditional



FS-WA00 / 10

Backing	PE foam, enhanced adhesion
Dimensions	45 × 32 mm
Gel	liquid
Packaging	50/pouch, 500/box, 3000/case


Premium Easitab, particularly for demanding conditions, 10 pcs./card, X-ray / MR conditional



RT34

Backing	Paper/PET
Dimensions	34 × 22 mm
Gel	solid
Packaging	100/pouch, 500/box, 5000/case

Easitab, 10 pcs./card, X-ray / MR conditional



FS-RB4 / 5

Backing	PE foam, enhanced adhesion
Dimensions	44 × 28 mm
Gel	liquid
Packaging	50/pouch, 1500/box, 6000/case

With adapter for 4 mm banana plug, 5 pcs./card



RT74

Backing	Paper/PET
Dimensions	34 × 26 mm
Gel	solid
Packaging	100/pouch, 500/box, 5000/case

Easitab, 10 pcs./card, X-ray / MR conditional

Monitoring

Adhere even better

Skintact® electrodes for monitoring ECG are available with standard adhesive or with reinforced adhesive – product names starting with „FS-“.

Easy to handle

An easy-to-open pouch and the peel-off aid for easy use.

Water repellent

All foam ECG electrodes are impermeable to liquids and thus protect the sensor.



FS-521

Backing PE foam, enhanced adhesion
Dimensions 50 × 35 mm
Gel solid
Packaging 30/pouch, 1500/box, 6000/case

Impermeable to liquids, ideal for short-term monitoring and for use in the operating room



CT-601

Backing transparent clear tape
Dimensions Ø 50 mm
Gel solid
Packaging 30/pouch, 1500/box, 6000/case

Skin reaction easily and quickly detectable



F-601

Backing PE foam
Dimensions Ø 50 mm
Gel solid
Packaging 30/pouch, 1500/box, 6000/case

Impermeable to liquids, ideal for short-term monitoring and for use in the operating room



T-60

Backing microporous tape
Dimensions Ø 60 mm
Gel liquid
Packaging 30/pouch, 1200/box, 4800/case

More comfort for patients due to breathable backing material, ideal for long-term monitoring



F-55

Backing PE foam
Dimensions Ø 55 mm
Gel liquid
Packaging 30/pouch, 1200/box, 4800/case

Impermeable to liquids, ideal for short-term monitoring and for use in the operating room, large adhesive area

Alternative: F-50 with 50 mm diameter



T-601

Backing microporous tape
Dimensions Ø 50 mm
Gel solid
Packaging 30/pouch, 1500/box, 6000/case

More comfort for patients due to breathable backing material, ideal for long-term monitoring



F-TB1

Backing PE foam
Dimensions 55 × 50 mm
Gel solid
Packaging 30/pouch, 1500/box, 6000/case

Impermeable to liquids, ideal for short-term monitoring and for use in the operating room



W-601

Backing soft cloth
Dimensions Ø 50 mm
Gel solid
Packaging 30/pouch, 1500/box, 6000/case

More comfort for patients due to breathable backing material, ideal for long-term monitoring

Universal

Lower costs

Our Skintact® ECG electrodes on cards of 6 or 10 pieces offer the usual Leonhard Lang quality at reasonable prices.

Fast and comfortable

These electrodes come on cards of 6 or 10 pieces. A perforation in the middle allows one card to be separated into two handy strips.

Feel free to leave the gloves on

The label material overlaps the electrode and creates a lift tab which allows the clinical staff to handle the electrode easily even when wearing surgical gloves.



FS-RG1 / 10

Backing PE foam, enhanced adhesion
Dimensions 41 × 32 mm
Gel solid
Packaging 50/pouch, 1500/box, 6000/case

Card of 10, also available on card of 6 (FS-RG1 / 6)



FS-RG / 6

Backing PE foam, enhanced adhesion
Dimensions 41 × 32 mm
Gel liquid
Packaging 30/pouch, 1200/box, 4800/case

Card of 6, also available with standard adhesion (F-RG / 6)



FS-TC1 / 10

Backing PE foam, enhanced adhesion
Dimensions 41 × 35 mm
Gel solid
Packaging 50/pouch, 1500/box, 6000/case


Card of 10, also available on card of 6 (FS-TC1 / 6)



FS-TF / 6

Backing PE foam, enhanced adhesion
Dimensions 51 × 43 mm
Gel liquid
Packaging 30/pouch, 1200/box, 4800/case

Card of 6, also available with standard adhesion (F-TF / 6)



MR-conditional & Radiotranslucent

■ Carbon snaps save time and costs

MR-conditional Skintact® ECG electrodes are equipped with carbon snaps. This saves time and costs: If the electrodes have already been applied, they do not have to be removed before X-ray or MRI examinations.

Skintact® Tab'n'Snap – two connection systems

With two integrated connection systems, Skintact® Tab'n'Snap electrodes ensure seamless patient transfer without the need for electrode replacement.



F-40C

Backing PE foam
Dimensions Ø 40 mm
Gel liquid
Packaging 30/pouch, 1200/box, 4800/case

X-ray / MR conditional for adults and pediatrics



FS-521C

Backing PE foam, enhanced adhesion
Dimensions 50 × 35 mm
Gel solid
Packaging 30/pouch, 1500/box, 6000/case

X-ray / MR conditional for adults



F-401C

Backing PE foam
Dimensions Ø 40 mm
Gel solid
Packaging 30/pouch, 1500/box, 6000/case

X-ray / MR conditional for adults and pediatrics



FS-RGC /10

Backing PE foam, enhanced adhesion
Dimensions 41 × 32 mm
Gel liquid
Packaging 50/pouch, 1500/box, 6000/case

X-ray / MR conditional for adults and pediatrics, card of 10



F-50C

Backing PE foam
Dimensions Ø 50 mm
Gel liquid
Packaging 30/pouch, 1200/box, 4800/case

X-ray / MR conditional for adults; available also with enhanced adhesion (FS-50C)



FS-WA10C

Backing PE foam, enhanced adhesion
Dimensions 45 × 32 mm
Gel liquid
Packaging 50/pouch, 500/box, 3000/case

Premium ECG electrode with Tab'n'Snap connection for adults, card of 10



F-601C

Backing PE foam
Dimensions Ø 50 mm
Gel solid
Packaging 30/pouch, 1500/box, 6000/case

X-ray / MR conditional for adults



W-130C /3

Backing soft cloth
Dimensions Ø 30 mm
Gel solid
Packaging 30/pouch, 1500/box, 6000/case

X-ray / MR conditional for all patient groups, adhesive area fully covered with soft hydrogel

Pediatrics & Neonates

Size and adhesive for pediatrics

Skintact® ECG electrodes for pediatrics and neonates are available in smaller sizes to meet the specific needs of smaller patients.

Offset connection


Pediatric offset electrodes, ensuring optimal positioning and minimizing the risk of artifacts – available with snap connection or a 4 mm banana adapter.

Pre-wired


For neonates, we offer pre-wired electrodes to minimize the risk of pressure necrosis.




© New Africa – stock.adobe.com




F-261
Backing PE foam
Dimensions Ø 26 mm
Gel solid
Packaging 30/pouch, 1500/box




F-VK01
Backing PE foam
Dimensions 45 × 30 mm
Gel liquid
Packaging 50/pouch, 500/box, 3000/case
Premium electrode with offset connection, available also with adapter for 4 mm banana plug (F-VK04)




F-301
Backing PE foam
Dimensions Ø 30 mm
Gel solid
Packaging 30/pouch, 1500/box




T-401
Backing microporous tape
Dimensions Ø 40 mm
Gel solid
Packaging 30/pouch, 1500/box
Breathable backing material




F-40
Backing PE foam
Dimensions Ø 40 mm
Gel liquid
Packaging 30/pouch, 1200/box



W-401
Backing soft cloth
Dimensions Ø 40 mm
Gel solid
Packaging 30/pouch, 1500/box
Breathable backing material




F-401
Backing PE foam
Dimensions Ø 40 mm
Gel solid
Packaging 30/pouch, 1500/box



W-130C /3
Backing soft cloth
Dimensions Ø 30 mm
Gel solid
Packaging 30/pouch, 1500/box, 6000/case
X-ray / MR conditional for all patient groups, soft hydrogel over the whole adhesive area fully covered with gel



PD50-SF (150)
Backing soft cloth
Dimensions 24 × 23 mm
Gel solid
Packaging 3/pouch, 150/box
With pre-attached 50 cm copper leads (red/yellow/black) with safety plug



PD50-F4C
Backing soft cloth
Dimensions 24 × 23 mm
Gel solid
Packaging 3/pouch, 150/box
With pre-attached 50 cm carbon leads (red/yellow/green) with connector for 4 mm banana plug; suitable for use during x-ray



PD50-SFC
Backing soft cloth
Dimensions 24 × 23 mm
Gel solid
Packaging 3/pouch, 150/box
With pre-attached 50 cm carbon leads (red/yellow/green) with safety plug; suitable for use during x-ray



PD10-F4
Backing soft cloth
Dimensions 24 × 23 mm
Gel solid
Packaging 10/pouch, 150/box
With pre-attached 10 cm copper leads with connector for 4 mm banana plug

Defibrillation Electrodes

When every second counts, patients around the world can rely on our Skintact® Defibrillation Electrodes.

In emergency cases not only the clinician's skills are required but also reliable ranges of life saving products like defibrillation electrodes. These products must perform dependably, even under adverse conditions as every second may count. Leonhard Lang's engineering team has developed reliable and safe multifunctional defibrillation electrodes. They can be used for external defibrillation, non-invasive pacing, synchronized cardioversion and monitoring.




Advantages of Skintact® Defibrillation Electrodes

- Unique overlapping gel technology, developed by Leonhard Lang, prevents tin edge exposure and helps to minimize the risk of burns.
- Secure adhesion even on perspiring skin.
- Maximum contact area combined with low impedance helps to reduce thermal and mechanical stress on the skin.
- Clear, two-color instructions for the correct placement of the electrodes.
- Multifunctional – one electrode for different applications.
- Radiotranslucent design minimizes visibility in X-ray imaging.
- Each set is subject to a comprehensive 100% fit and continuity inspection.

Information

All our Skintact® products are tested according to applicable standards, ensuring maximum safety and reliable performance.

Skintact® Easibeat Multifunction Defibrillation Electrodes meet the relevant clauses of IEC 60601-2-4 for Cardiac Defibrillators, and do comply with ISO 10993 for biocompatibility and the European Medical Device Regulation (MDR).



Multifunctional Radiotranslucent Defibrillation Electrodes for Adults

Used for external defibrillation, non-invasive pacing, synchronized cardioversion and monitoring for adults and children.


DF20N

Compatible with **Physio – Control:** Lifepak 12, Lifepak 15, Lifepak 20/20e, Lifepak 500, Lifepak 1000

Comen: S1, S5, S8

Cable length 1,2 m

Color code ■




DF20NC

Compatible with **Physio – Control – preconnect:** Lifepak 12, Lifepak 15, Lifepak 20/20e, Lifepak 500, Lifepak 1000

Cable length 1,2 m

Color code ■




DF27N

Compatible with **Philips Heartstart:** FR2, FR2+, FR3 (electrode self-test is not supported) With adapter cable M3508A or M3507A and adapter 05-10200: Philips Heartstart MRx, XL, XL+, Efficia DFM100, Intrepid

Cable length 2,0 m

Color code ■




DF27NC

Compatible with **Philips Heartstart – preconnect:** FR2, FR2+, FR3 (electrode self-test is not supported) With adapter cable M3508A or M3507A and adapter 05-10200: Philips Heartstart MRx, XL, XL+, Intrepid

Cable length 2,0 m

Color code ■




DF28N

Compatible with **Zoll:** Series M, E, R and X With adapter cable 8000-0308-01/02, 8009-0749, 8300-0783, 8300-000676

Cable length 1,2 m

Color code ■




DF28NC

Compatible with **Zoll – preconnect:** Series M, E, R and X With adapter cable 8000-0308-01/02, 8009-0749, 8300-0783, 8300-000676

Cable length 1,2 m

Color code ■



DF26N

Compatible with **Philips:** With adapter cable M3507A: Philips Heartstart MRx, XL, XL+, Efficia DFM100, Intrepid

Cable length 1,2 m



DF43N

Compatible with **Nihon Kohden:** With pad adapter JC-755V: TEC-5500 Series, TEC-7600 Series, TEC-7700 Series With pad adapter JC-855V: TEC-8300 Series

Cable length 2,0 m



Packaging
All Skintact® Easibeat Multifunction Defibrillation Electrodes
1/pouch 10/box 60/case

Multifunctional Defibrillation Electrodes for Pediatrics

Used for external defibrillation, non-invasive pacing, synchronized cardioversion and monitoring of pediatric patients.

DF31G

Compatible with **Physio – Control:** Lifepak 12, Lifepak 15, Lifepak 20/20e

Comen: S1, S5, S8


Cable length 1,2 m



DF30G

Compatible with **Zoll:** Series M, E, R and X With adapter cable 8000-0308-01/02, 8009-0749, 8300-0783, 8300-000676

Cable length 1,2 m



DF32G

Compatible with **Philips Heartstart:** With adapter cable M3508A or M3507A and adapter 05-10200: Philips Heartstart MRx, XL, XL+, Efficia DFM100, Intrepid


Cable length 1,2 m



DF37G

Compatible with **Nihon Kohden:** With electrode adapter JC-755V: TEC-5500 Series, TEC-7600 Series, TEC-7700 Series With electrode adapter JC-855V: TEC-8300 Series

Cable length 2,0 m




Neutral Electrodes


Especially in electrosurgical procedures, it is important that electrodes and accessories provide absolute safety.

Maximising the full potential of advanced electrosurgery in today's market, not only requires the surgeon's skills but also a reliable range of active and passive accessories. These accessories should offer safety and comfort without compromise to both patient and clinician. Leonhard Lang's range of Skintact® Cool Contact electrosurgical neutral electrodes and pencils provide the highest degree of safety.



Advantages of Skintact® Neutral Electrodes

- Skintact Cool Contact neutral electrodes are available for adults and pediatric patients, with or without a preattached cable.
- Developed by Leonhard Lang, the unique overlapping gel technology eliminates aluminum edge exposure and helps to minimize the risk of burns.
- The rounded leading-edge geometry of the universal neutral electrodes enables more efficient current distribution, helping to maintain a cooler interface between the patient and the neutral electrode.
- The peel-off aid allows easy removal from the liner or the patient's skin.
- The flexible, thin foam backing of the RS and RO series as well as the low-profile fabric backing of the WR series are impermeable to liquids. They conform to the body's contours and provide reliable protection and safety.
- Duplex labels reduce the documentation effort for clinical staff.




Information

All our Skintact® products are tested according to applicable standards, ensuring maximum safety and reliable performance.

Skintact® Cool Contact Electrosurgical Neutral Electrodes comply with the applicable clauses of the IEC 60601-2-2 standard for high-frequency surgical devices, the ISO 10993 standard for biocompatibility and the European Medical Device Regulation (MDR).

Neutral Electrodes



RO21
RO21/5 *

Size universal
Weight limit > 5 kg
Conductive area 107 cm²




WR21
WR21/5 **

Size universal
Weight limit > 5 kg
Conductive area 90 cm²



RS25
RS25/5 *

Size adult
Weight limit > 15 kg
Conductive area 118 cm²



RS27
RS27/5 *

Size adult
Weight limit > 15 kg
Conductive area 120 cm²



RS26

Size pediatric
Weight limit 5–15 kg
Conductive area 66 cm²
Power maximum 120 watts

Neutral Electrodes with cables



RO21A30

Size universal
Weight limit > 5 kg
Conductive area 107 cm²
Cable length 3 m



WR21A30

Size universal
Weight limit > 5 kg
Conductive area 90 cm²
Cable length 3 m
upon request we also offer 5 m



RS25A30

Size adult
Weight limit > 15 kg
Conductive area 118 cm²
Cable length 3 m
upon request we also offer 5 m



RS27A30

Size adult
Weight limit > 15 kg
Conductive area 120 cm²
Cable length 3 m



RS26A30

Size pediatric
Weight limit 5–15 kg
Conductive area 66 cm²
Power maximum 120 watts



RSW23A30 ***

Size newborn
Weight limit 0,45–4,99 kg
Conductive area 43 cm²
Current maximum 350 mA

Pencils & Tip Cleaner

- Single use eliminates risk of cross infection.
- Sterile, easy to open packaging.
- More cost effective than cleaning reusable versions.



HCP-01

Sterile Fingerswitch Pencil with standard blade

Cable length 3 m
Connector 3 Pins at 4 mm each
Packaging 50/box



HCP-02

Sterile Fingerswitch Pencil with standard blade


Cable length 3 m
Connector 1 Pin at 5 mm
Packaging 50/box



HCP-05

Sterile Fingerswitch Pencil with standard blade

Cable length 5 m
Connector 3 Pins at 4 mm each
Packaging 50/box



TC 1001

Sterile Tip Cleaner

Dimensions 50 × 50 mm
Packaging 100/box

Reusable Cables

Reusable neutral electrode cable with embossed lot number for easy tracking.



NEK300A
NEK500A

Cable length 3 m
5 m
Connector REM™ (Valleylab™)
Packaging 1/pouch



NEK300D
NEK500D

Cable length 3 m
5 m
Connector 6,3 mm jack plug (Standard)
Packaging 1/pouch



NEK300H
NEK500H

Cable length 3 m
5 m
Connector 6,3 mm jack plug (Erbe)
Packaging 1/pouch



Packaging
Skintact® Cool Contact
Neutral Electrodes

1/pouch 50/box
* 5/pouch 50/box
** 5/pouch 100/box
*** 1/pouch 10/box

