



Patient **Support** Systems
MEDICAL • DESIGN • INNOVATION

The Alliana Seat Cushion

Affordable pressure care



- ✓ Soft PVC Air Cells
- ✓ Adjustable comfort
- ✓ All seams fully welded
- ✓ Belgian two way stretch covers
- ✓ Zip closure for easy access & cleaning

Specifications

Alliana Seat Cushion

Dimensions Standard	L 48 x W 44 x H 6cm
Cell Type	56 x PVC Soft Single Mould Cells
Top Cover	Dahlia Bi-Elastic 2 way stretch - welded
Base Standard	Dahlia Bi-Elastic 2 way stretch - welded

Order Code:

ALLI-STC-44x48x6

Inflation	Hand pump and valve	
Infection Control	All seams fully welded - waterfall over zip	
Base Attachment Standard	Rear Velcro Straps	
Warranty	2 years Soft Goods	ARTG: 175810


Unit C3, 15 Narabang Way Belrose, Sydney NSW 2085
t: 02 9844 5456 f: 02 9844 5445 e: info@patientsupportsystems.com

www.patientsupportsystems.com

The Alliana Air Alternating Cushion


BlueTooth 3 in 1 Therapy Control

- ✓ Six soft PU air cells
- ✓ Alternation, pulsation and static modes
- ✓ Precise bluetooth control from phone
- ✓ 14 hour battery
- ✓ Zip closure for easy cleaning



Alternation A

Alternation B



Pressure Mapping 6/5/21

Cerebral Palsy Alliance

Male 92kg

Alliana Alternating BlueTooth Cushion

(Wheelchair not included)

Dimensions: L 45x W 44 x H 7.5cm

Warranty: 2 years

ARTG: 175810

Call 1800 432 247 to order yours today