

Viable Surface Sampling Worksheet for Category 3 CSP in ISO Class 5 PEC

Sampling Date: _____

Type of Media	Manufacturer	Lot Number	Expiration Date
TSA 2/ Lecithin & Polysorbate 80			
SDA Sabouraud Dextrose Agar w/ Lecitin & Polysorbate 80			

TSA (30°-35° Incubator ID: A)	SDA (20°-25° Incubator ID: B)
Date: _____ & Time: _____ Initials of Staff placing in Incuabor: _____	Date: _____ & Time: _____ Initials of Staff placing in Incuabor: _____

Lot Number	# Observed CFU								Action Levels	PASS Below Action Level	FAIL Equal or Above Action Level	Comments about specific plate made below
	TSA 30°-35° for 3 days			SDA 20°-25° for 5 days								
	1	2	3	1	2	3	4	5				
									>3			
									>3			
									>3			
									>3			
									>3			
									>3			
									>3			
									>3			
									>3			
									>3			
									>3			
									>3			
									>3			
									>3			
									>3			
									>3			
									>3			
Date									© 2006-2026 Bayview Pharmacy Inc. All rights reserved.			
Initials of Staff Reading												