

Ordering

500885 - 500 g

Applications

Suspensions

Product Highlights

SyrSpend® SF Alka dry is a pre-measured alkaline suspending base for reconstitution to compound elegant oral liquid dosage forms with a pleasant taste. It comes in a validated, graduated dispensing container, allowing for an easy three-step compounding process directly into the dispensing container. This unique vehicle uses an active suspending technology to offer more accurate dosing, by stabilising the suspension at rest, and actively preventing 'caking' of insoluble ingredients at the bottom of the bottle.

- Ready-to-use, all-in-one suspending and sweetening liquid vehicle
- Compatible with an extensive range of APIs
- Buffered to pH >7.0 for maximum stability of acid-labile APIs
- Low osmolality (<50 mOsmol/kg)
- Preservative-free

Ingredients

Calcium carbonate, Food starch modified, Sucralose

Specifications

SyrSpend SF ALKA Dry is a white, free-flowing powder

- pH: >7.0 (when reconstituted)

Formulas

#FA21488 Esomeprazole 3 mg/mL Oral Suspension in SyrSpend SF ALKA Dry

#FA21312 Furosemide 10 mg/mL Oral Suspension in SyrSpend SF ALKA Dry

#FA21389 Pantoprazole 2 mg/mL Oral Suspension in SyrSpend SF ALKA Dry

Formulas are available for download for FACTS members: fagronacademy.us

Studies

Fagron scientific literature fagronacademy.us/facts/literature