




Mississippi WIC Enrollment Codes & Detailed Enrollment Overview



Enrollment Codes

Client Name: Mississippi WIC

What is a Pacify enrollment code?

This is a unique code enrollees enter during the Pacify app sign up process for free access to Pacify. This code allows Pacify to link that member to the Mississippi WIC services. This information will be included in the weekly enrollment report provided to the Mississippi WIC team.

Need a new Pacify enrollment code created? Please contact support@pacify.com. If the code is for a new staff member please include the name, team and title of the staff member. If the code is for a new enrollment avenue please state the channel (e.g. newsletter). This usually takes 2 days.

Code	Code Description
MS011	District 1
MS012	District 2
MS013	District 3
MS024	District 4
MS025	District 5
MS026	District 6
MS037	District 7
MS038	District 8
MS039	District 9
MSWICAdmin	Staff

Code	Code Description
PCDBE1	Peer Counselor District 1
PCDBE2	Peer Counselor District 2
PCDBE3	Peer Counselor District 3
PCDBE4	Peer Counselor District 4
PCDBE5	Peer Counselor District 5
PCDBE6	Peer Counselor District 6
PCDBE7	Peer Counselor District 7
PCDBE8	Peer Counselor District 8
PCDBE9	Peer Counselor District 9



Detailed Enrollment Overview

STEP 1: Introduce Pacify

Pacify is a mobile app that provides video-enabled support from maternal health specialists - all at the touch of a button!



On-demand



Unlimited
calls



No cost



Bilingual

STEP 2: Assist in Account Creation

1. Help the participant download Pacify. It's available for free on both Apple and Android devices. Then, have them open Pacify.
2. Help the participant create an account. Select "Create an account" or "Don't have an account" at the bottom of the login screen. Then, have them enter an email address and password.
3. On the next screen, have the participant enter the appropriate sign-up code.
4. Assist the participant while they complete their profile.

STEP 3: Encourage a Welcome Call

1. Explain that lactation calls are video-enabled, just like FaceTime.
2. Instruct the participant to start a "test" call, by pressing the Lactation Consultant button.
3. Have the participant initiate the video call, ask any questions, and then leave feedback when done.

