


BEFORE & AFTER⁺




Female 30y • 1 session of 3mm H/P 100W, 60000J • Immediately After



Female 20y • 1 session of 3mm H/P 100W, 60000J • Immediately After



Male 20y • 2 session of SHEP / 60KJ / 100WATT • 2 month



Female 30y • 1 session of 7mm H/P 100W, 60000J • Immediately After



Female 30y • 1 session of 7mm H/P 90W, 60000J • Immediately After

KEY BENEFITS⁺



FAT REDUCTION

Precisely targets and heats fat cells beneath the skin to reduce stubborn fat in targeted areas.



BODY CONTOURING

Enhances your natural curves and creates a more sculpted, balanced silhouette.



SKIN TIGHTENING

Stimulates collagen production to improve skin elasticity and firmness for a smoother, youthful look.



FACIAL LIFTING

Helps define facial contours, lift sagging areas, and improve overall skin firmness.



CELLULITE IMPROVEMENT

Reduces the appearance of cellulite by addressing its root cause and smoothing the skin surface.

WHY CHOOSE MIGLOU^{TM+}

- ✓ Non-invasive, no surgery
- ✓ No downtime, return to daily life immediately
- ✓ Precisely targets fat without damaging surrounding tissues
- ✓ Customizable treatments for your unique body
- ✓ Safe, comfortable, and clinically designed
- ✓ Visible results with continued improvement over time

Let Me Glow,
MIGLOUTM



Signature
HIRONIC Ambassador
YANO SHIHO

FOCUSED ENERGY
DESIGNED GLOW

Let Me Glow,
miglou™

PRECISION-DRIVEN SKIN TIGHTENING



01 2.45 GHz UHF WITH POWER COOLING

This technology efficiently delivers energy to the desired depth, effectively transferring heat to the targeted area while protecting the skin through cooling while efficiently transmitting energy.

02 SUPPORTS SKIN CONTOURING AND TIGHTENING


Targets fat cells and stimulates collagen production, leading to fat reduction and collagen regeneration.

03 SMART HANDPIECES BY PENETRATION DEPTHS

MIGLOU™ offers customized treatment for each patient and enhanced convenience of operation. Choose our various handpieces for different indications for use.

FIM* SMART MATCHING⁺

- ✓ Focused Intensive Microwave technology enhances the concentration and efficiency of delivered energy
- ✓ MIGLOU's proprietary, patented technology precisely delivers energy to targeted focal zones with high accuracy



CONTACT & MOVING SMART SENSOR⁺



SAFER TREATMENT WITH SMART SENSOR TECHNOLOGY

- Built-in handpiece safety system helps prevent overheating
- Detects both contact and movement for controlled use
- Contact cooling system helps protect the skin and provide comfort treatment


TREATMENT AREAS⁺

FACIAL TREATMENT

- Double Chin Reduction
- Skin Tightening
- Facial Lifting

BODY TREATMENT

- Abdomen
- Love Handles
- Thighs
- Upper Arms
- Calves
- Back



FREQUENTLY ASKED QUESTIONS⁺

Q1. How long does the treatment take?

A. Treatment takes about 15 minutes, depending on the treatment area.

Q2. Is it painful?

A. It is designed for a comfortable treatment with no numbing cream required.

Q3. Is there downtime?

A. You can generally return to normal daily activities.

Q4. What areas can be treated?

A. Face, abdomen, thighs, love handles, calves, upper arms, back.

Q5. How many sessions?

A. Face: 3–4 sessions every 3–4 weeks.
Body: 5–6 sessions every 4 weeks.

Q6. What results can you expect?

A. Improved contouring, tightening, fat reduction, and smoother skin texture. Results continue to develop as collagen production increases over time.