

ULTRAFLO-C

Capillary Blood Collection System

The ideal combination of controlled blood collection and sample container to avoid sample loss or contamination

Exact measurement and stabilization of the required amount of blood via integrated end-to-end plastic capillary coated with anti-coagulant.

Reduction of the infection risk thanks to unbreakable plastic capillaries

Clear marking of expiry date and production batch on each sample container



Capillary
volume
100µl
200µl

Preparation
EDTA-K2
EDTA-K2

Cap color
lavender
lavender

REF No.
681005
682005

Compare to
ASP(RAM)
07 6011
07 7051

MADE IN GERMANY FOR



INNOVATIVE MED TECH

866-560-1820
support@InnovativeMedTech.com
www.InnovativeMedTech.com