

NitriDerm[®] 100

Nitrile Exam Gloves Series 164

Engineered for exceptional value, this nitrile medical glove combines superior comfort and performance. Crafted from a soft, low-modulus synthetic material, it provides excellent elasticity and is easy to don. Conveniently packaged in a standard 100-count box.

Characteristics: Non-Sterile, Powder-Free, Blue Color



SIZE	ITEM NO.
XS	164050
S	164100
M	164200
L	164300
XL	164350

XS-XL : 100 Gloves/Box, 10 Boxes/Case



Tested for use with
Chemotherapy Drugs & Fentanyl Citrate

Features:

- Non-Latex = No Risk of Latex Allergens
- Textured Finish = Improved Wet/Dry Grip
- Low Modulus = a Softer, More Comfortable Fit

Manufactured in compliance with multiple international standards, including ASTM D6319, ASTM D5151, ASTM F1671, and ASTM D6978

Tested for use with Chemotherapy Drugs as per ASTM D6978

Standard Practice for Assessment of Medical Gloves to Permeation by Chemotherapy Drugs

Minimum Breakthrough Detection Time in Minutes:

Carmustine (3.3 mg/ml)	14.1	Methotrexate (25.0 mg/ml)	>240
Cisplatin (1.0 mg/ml)	>240	Mitomycin C (0.5 mg/ml)	>240
Cyclophosphamide (20.0 mg/ml)	>240	Oxaliplatin (5 mg/ml)	>240
Dacarbazine (10 mg/ml)	>240	Paclitaxel (6.0 mg/ml)	>240
Doxorubicin Hydrochloride (2.0 mg/ml)	>240	Thiotepa (10.0 mg/ml)	57.2
Etoposide (20.0 mg/ml)	>240	Vincristine Sulfate (1.0 mg/ml)	>240
Fluorouracil (50.0 mg/ml)	>240	Fentanyl Citrate (100 mcg/2ml)	>240

Please note that the following drugs have low permeation times: Carmustine 14.1 minutes, Thiotepa 57.2 minutes.

Caution: In compliance with ASTM D6978, the testing conditions used are intended to approximate the worst-case conditions for clinical use. Testing was conducted on single layer glove material. It is the users' responsibility to determine the applicability of these gloves for their intended use with chemotherapy drugs.

Tensile Strength & Elongation	Before Aging	After Accelerated Aging
Tensile Strength (Mpa)	25	29
ASTM Requirement Min. (Mpa)	14	14
Elongation (%)	530	460
ASTM Requirement Min. (%)	500	400

Average Length	9.5 in ♦ 240 mm
Average Palm Thickness	2.4 mil ♦ 0.06 mm
Average Finger Thickness	2.8 mil ♦ 0.07 mm



Innovative Healthcare Corporation is certified to ISO 13485:2016 QMS for medical devices.



INNOVATIVE
HEALTHCARE
CORPORATION

www.ihc-solutions.com

1-800-272-1533

View this product
on our website:

