

RX1-Mini

Mobile Cardiac Telemetry (MCT) & Events



New

The RX-1 Mini provides Multichannel MCT and Event Monitoring in a single, compact device—for complete and responsive patient care.



Built-in 4G. No special hardware or cellphone needed.

Faster, more comprehensive clinical reports that make a difference for your patients.

- Urgent notifications & episode reports
- Daily Summary & End of Service reports
- Heart rate histogram on MCT report
- A-Fib onset and offset on MCT summary report
- AF burden and rate charts
- P-wave enhancement available on portal
- Full disclosure

Built with the patient in mind, the RX-1 Mini is easy to wear and use without compromising on accurate diagnostic data.

Onboard cellular means no more pairing with other devices and simpler day-to-day management for patients being monitored in long-term cardiac studies.

One-button operation and long battery life make it easy for the wearer. The predetermined algorithms and 2 channels of detection ensure clean and accurate data.

Clinical benefits

Multichannel

Get it right the first time. All reports contain 2-channels of ECG data, ensuring accurate and comprehensive reporting.

Quicker turnarounds

Using advanced algorithms, arrhythmias are identified faster. Signal noise is greatly reduced, leading to faster analysis.

Long battery life

Up to 5 days of device battery life on a single 1-hour charge.

Universal lead connectors

No need for expensive, proprietary electrodes. Deliver cost-effective care and provide electrodes tailored to each patient's skin type.

Patient benefits

Lightweight & comfortable

At 1.6 oz, the RX-1 Mini is small and unobtrusive, meaning patients can continue their daily routines uninterrupted.

Simple one-touch operation

No complicated interfaces required—once the device is recording just hold the button for 5 seconds to mark an event.

Durable & showerproof

Built to handle everyday drops and spills. Patients can go about their daily routines without compromising data integrity or study results.

Wireless data transfer

No special hardware or cellphone needed. Data is automatically sent to the cloud using built-in 4G.

Contact us

Cardiac Monitoring Service
1300 Quail Street, Suite 209 Newport Beach, CA 92660

T: (800)-267-5248
E: info@cmscardioscan.com

