

# Studio Bariatric

## 4 4 World's first cloud connected therapy surface



Precise control and monitoring  
iPhone, iPad or Windows PC

## 4 Hyperstretch coverlet standard with sacral cell stabilisation

## 4 Electronic valves for 25% better therapy and proven wound healing

## 4 Remote alerts for total peace of mind

Static inflatable side cells for additional  
transfer support



Independent Pressure Mapping by Cerebral Palsy Alliance May 2021  
Male subject 92kg sacrum



**STUDIO CARE**

# TheraCloud

## Bariatric Specifications



### Control Unit:

|                                                     |                      |
|-----------------------------------------------------|----------------------|
| Dimensions:                                         | W 27 x H 21 x D 10cm |
| Weight                                              | 2.8kg                |
| Mode of Operation                                   | Non-Continuous       |
| Power Rating                                        | 240V 50Hz 12VA       |
| Transport Function                                  | Hose End Cap         |
| Air Flow                                            | 18 Litres per minute |
| Auto Startup                                        | 3                    |
| Auto Sensor for Patient Weight Profile and Position | 3                    |
| Comfort Override Control                            | 3                    |
| Self Diagnostics                                    | 3                    |
| Audible/Visual Alarm with email alerts              | 3                    |
| Mute Function                                       | 3                    |
| Static Mode                                         | Timeout - 60 minutes |
| Care Mode                                           | Timeout - 20 minutes |
| Incline mode with Auto Tilt Sensor                  | 3                    |
| Auto Anti Tamper - 2 Second Unlock Function         | 3                    |
| CE Certification and C-Tick Listed N27882           | ✓                    |

### Warranty

2 Years Pump  
2 Year Soft Goods

ARTG: 175810

### Mattress Replacement

Dimensions L 200 x W 120 x H 23cm Alternating Cell Type  
Safety Cell in Cell Overall Cells Count 24 - (3 static  
Head/Side Cells included) Quick release valves for last foot  
zone cells 6 Cell Material Pure TPU Top Cover Sioen 2 way  
stretch - welded Base/s Nylon TPU welded Cell Cycle 1 in 2  
Alternating Mode 10 minutes per alternation Infection  
Control Welded top cover, handles and base Positioning  
Handles 2 aside Base Attachment Non Slip Base Material

### Hyperstretch Coverlet

Dimensions L 200 x W 120 x H 1cm Attachment Base  
Surround Zip Cover Material High Performance Micromesh  
Foam Dunlop Marathon Hypersoft HS18-35 with Ultrafresh



STUDIO CARE