

# Certificate of Analysis

## LIM Kinase 1, active

(Recombinant enzyme expressed in Sf21 insect cells) Item # 14-656, 14-656-K, 14-656M Parent Lot # WAB0453

The data presented in this document apply to the parent lot shown above and to all pack sizes derived from subsequent vialling runs of this parent lot. An alphabetical suffix after the parent lot number is used to denote each vialling run

**Product Description:** *N*-terminal 6Histagged, recombinant, human LIM Kinase 1, amino acids 285–638, expressed by baculovirus in Sf21 insect cells. Purified using Ni<sup>2+</sup>-NTA agarose. Activated using ROCKII (cat# 14-451) and repurified using Q Sepharose. Purity 70% by SDS-PAGE and Coomassie blue staining. MW = 41.2kDa.

Specific Activity (Parent lot# WAB0453): 1813U/mg, where one unit of LIM Kinase 1 activity is defined as 1nmol phosphate incorporated into 0.634mg/ml cofilin 1 per minute at  $30^{\circ}C$  with a final ATP concentration of  $100\mu M$ .

**Formulation: 0.224mg/ml** of enzyme in 50mM Tris/HCl pH7.5, 300mM NaCl, 0.1mM EGTA, 0.03% Brij-35, 270mM sucrose, 1mM benzamidine, 0.2mM PMSF, 0.1% 2-mercaptoethanol. Frozen solution.

**Storage and Stability:** On receipt of material store at -70°C. Unopened reagent is stable for a minimum of 1 year from date of shipment when stored at recommended storage temperature. Avoid repeat freeze/thaw cycles. For maximum recovery of product, centrifuge original vial prior to removing the cap.

**Handling Recommendations:** Rapidly thaw the vial under cold water and immediately place on ice. Aliquot unused material into pre-chilled micro-centrifuge tubes and immediately snap-freeze the vials in liquid nitrogen prior to re-storage at -70°C.

# FOR IN VITRO RESEARCH USE ONLY NOT FOR USE IN HUMANS OR ANIMALS

97kDa -66kDa -

45kDa -


30kDa -

20.1kDa -


14.4kDa -

### **Quality Control Testing**

<u>Kinase Assay</u>: 0.75–16ng of this lot of enzyme phosphorylated 0.634mg/ml Cofilin 1 in the assay described on page two. Assay background was subtracted from the actual counts to yield the results shown below.



MS Tryptic Fingerprint: Confirmed product identity as LIM Kinase 1 with the translated native sequence listed on page three.





# Certificate of Analysis

### **Kinase Assay Protocol**

### Stock Solutions:

- 5 x Reaction Buffer: 40mM MOPS-NaOH pH7.0, 1mM EDTA.
- Cofilin 1: Use at a final concentration of 0.634mg/ml. Make a 3.17mg/ml stock. Add 5μl of stock per assay point.
- 3. LIM Kinase 1, active: Dilute with 20mM MOPS-NaOH pH7.0, 1mM EDTA, 0.01% Brij-35, 5% glycerol, 0.1% 2-mercaptoethanol, 1mg/ml BSA. Use 0.75–16ng per assay point.
- **4.** [ $\gamma$ -<sup>33</sup>P]ATP: 2.5 x magnesium acetate/[ $\gamma$ -<sup>33</sup>P]ATP cocktail: 25mM MgAc and 0.25mM ATP to which is added [ $\gamma$ -<sup>33</sup>P]ATP (specific activity approximately 500 800cpm/pmol as required.)

#### Assay Procedure (96 well plate format):

- 1. Add 5µl of 5 x reaction buffer per assay to wells.
- 2. Add 5µl of 3.17mg/ml cofilin 1.
- 3. Add 2.5µl (0.75–16ng) LIM Kinase 1, active.
- 4. Add 2.5µl of dH<sub>2</sub>O.
- 5. Add 10 $\mu$ l of diluted [ $\gamma$ -<sup>33</sup>P]ATP mixture.
- 6. Incubate for 15 minutes at 30°C.
- 7. Stop the reaction by adding 5µl of 3% phosphoric acid.
- 8. Transfer a 10µl aliquot onto the appropriate area of a **P30 Filtermat**.
- 9. Wash the filtermat three times for 5 minutes with 75mM phosphoric acid.
- 10. Wash the filtermat once for 2 minutes with methanol.
- 11. Transfer the filtermat to a sealable plastic bag and add 4ml of scintillation cocktail.
- 12. Read in a scintillation counter. Compare cpm of enzyme samples with cpm of control samples that contain all assay components plus 1µl of 30% phosphoric acid.

# **Certificate of Analysis**

#### **LIM Kinase 1 Sequence Information**

Protein Human LIM Kinase 1

**Tags** *N*-terminal 6His

Native sequence G10 of the recombinant protein is equivalent to G285 of human LIM Kinase 1

Accession number GenBank NM\_002314

#### Recombinant LIM Kinase 1 amino acid sequence:

```
1 MHHHHHHEFG SSARQKPVLR SCSIDRSPGA GSLGSPASQR KDLGRSESLR VVCRPHRIFR
61 PSDLIHGEVL GKGCFGQAIK VTHRETGEVM VMKELIRFDE ETQRTFLKEV KVMRCLEHPN
121 VLKFIGVLYK DKRLNFITEY IKGGTLRGII KSMDSQYPWS QRVSFAKDIA SGMAYLHSMN
181 IIHRDLNSHN CLVRENKNVV VADFGLARLM VDEKTQPEGL RSLKKPDRKK RYTVVGNPYW
241 MAPEMINGRS YDEKVDVFSF GIVLCEIIGR VNADPDYLPR TMDFGLNVRG FLDRYCPPNC
301 PPSFFPITVR CCDLDPEKRP SFVKLEHWLE TLRMHLAGHL PLGPQLEQLD RGFWETYRRG
361 ESG
```

#### Recombinant LIM Kinase 1 nucleotide sequence:

```
1 atgcatcatc accatcacca tgaattcggc agctctgccc ggcagaaacc tgtcttgagg
  61 agctgcagca tcgacaggtc tccgggcgct ggctcactgg gctccccggc ctcccagcgc
 121 aaggacctgg gtcgctctga gtccctccgc gtagtctgcc ggccacaccg catcttccgg
 181 ccgtcggacc tcatccacgg ggaggtgctg ggcaagggct gcttcggcca ggctatcaag
 241 gtgacacacc gtgagacagg tgaggtgatg gtgatgaagg agctgatccg gttcgacgag
 301 gagacccaga ggacgttcct caaggaggtg aaggtcatgc gatgcctgga acaccccaac
 361 gtgctcaagt tcatcggggt gctctacaag gacaagaggc tcaacttcat cactgagtac
 421 atcaagggcg gcacgctccg gggcatcatc aagagcatgg acagccagta cccatggagc
 481 cagagagtga gctttgccaa ggacatcgca tcagggatgg cctacctcca ctccatgaac
 541 atcatccacc gagacctcaa ctcccacaac tgcctggtcc gcgagaacaa gaatgtggtg
 601 gtggctgact tcgggctggc gcgtctcatg gtggacgaga agactcagcc tgagggcctg
 661 cggagcctca agaagccaga ccgcaagaag cgctacaccg tggtgggcaa cccctactgg
 721 atggcacctg agatgatcaa cggccgcagc tatgatgaga aggtggatgt gttctccttt
 781 gggatcgtcc tgtgcgagat catcgggcgg gtgaacgcag accctgacta cctgccccgc
841 accatggact ttggcctcaa cgtgcgagga ttcctggacc gctactgccc cccaaactgc
 901 cccccgagct tcttccccat caccgtgcgc tgttgcgatc tggaccccga gaagaggcca
961 tcctttgtga agctggaaca ctggctggag accctccgca tgcacctggc cggccacctg
1021 ccactgggcc cacagctgga gcagctggac agaggtttct gggagaccta ccggcgcggc
1081 gagagcggat ga
```

Reviewed and approved by site quality representative.

Unless otherwise stated in our catalogue or other company documentation accompanying the product(s), our products are intended for research use only and are not to be used for any other purpose, which includes but is not limited to, unauthorized commercial uses, in vitro diagnostic uses, ex vivo or in vivo therapeutic uses or any type of consumption or application to humans or animals.

© 2014 Eurofins Pharma Discovery Services UK Limited is an independent member of Eurofins Discovery Services