

Certificate of Analysis

Parkin, active

(Recombinant E3 ligase expressed in Sf21 insect cells) Item # 23-048, 23-048-K, 23-048M
Parent Lot # D11PP049N

The data presented in this document apply to the parent lot shown above and to all pack sizes derived from subsequent vialling runs of this parent lot. An alphabetical suffix after the parent lot number is used to denote each vialling run.

Product Description: *N*-terminal 6His-tagged, recombinant human Parkin full length, expressed by baculovirus in Sf21 insect cells. Purified using immobilized metal affinity chromatography.

Purity 85% by SDS-PAGE and Coomassie blue staining. MW = 55 kDa.

Activity (Parent lot# D11PP049N): This lot of Parkin is active and meets product specifications.

Formulation: 1.267mg/ml of enzyme in 50mM Tris/HCl pH7.5, 150mM NaCl, 0.03% Brij-35, 270mM sucrose, 1mM benzamidine, 0.2mM PMSF, 0.1% 2-mercaptoethanol. Frozen solution.

Storage and Stability: On receipt of material store at -70°C. Unopened reagent is stable for a minimum of 1 year from date of shipment when stored at recommended storage temperature. Avoid repeat freeze/thaw cycles. For maximum recovery of product, centrifuge original vial prior to removing the cap.


Handling Recommendations: Rapidly thaw the vial under cold water and immediately place on ice. Aliquot unused material into pre-chilled micro-centrifuge tubes and immediately snap-freeze the vials in liquid nitrogen prior to re-storage at -70°C.


FOR IN VITRO RESEARCH USE ONLY NOT FOR USE IN HUMANS OR ANIMALS

Quality Control Testing

<u>Assay</u>: This enzyme was titrated in a ubiquitination assay and the results normalised against the maximum signal.

<u>Protein Identity:</u> Confirmed identity as Parkin by mass spectrometry.





SDS-PAGE and Coomassie
Stain: Purity was assessed by
SDS-PAGE and Coomassie
blue staining using 3µg of
Parkin, active.



Certificate of Analysis

E3 Assay Protocol

Reagents:

- 1. UBE1, active (Item # 23-021)
- 2. UbcH7, active (Item # 23-047)
- 3. Parkin, active (Item # 23-048)
- 4. p38/JTV-1 (Item # 23-049)

- 5. 1x Reaction Buffer
- 6. Biotinylated-Ubiquitin
- 7. Stop Solution

Assay Outline:

All enzymes and reagents are diluted in the 1x reaction buffer (25mM MOPS pH 7.5, 0.01% Tween 20, 5mM $MgCl_2$).

Parkin is incubated with 25mM MOPS pH 7.5, 0.01% Tween 20, 5mM MgCl $_2$, 10 μ M ATP, 10nM UBE1, 500nM UbcH7, 50nM p38/JTV-1, and 2 μ M biotinylated-ubiquitin. The reaction is initiated with the addition of biotinylated-ubiquitin. After 30 minutes at room temperature the reaction is terminated by the addition of 25mM MOPS pH7.5 containing 125mM EDTA, 150mM NaCl, and 0.05% Tween 20. Reaction products are separated by capture onto a microplate coated with anti-c-Myc antibody and washing with PBS containing 0.05% Tween-20. Parkin activity is measured by detection of bound ubiquitin via electrochemiluminescence.



Certificate of Analysis

Parkin Information

Protein human Parkin

<u>Accession number</u> GenBank NM_004562

Alternative Names Parkinson juvenile disease protein 2, PARK2, PRKN

Key Facts Parkin is an E3-ubiquitin ligase belonging to the RBR (RING – In Between RING –

RING) family. Mutations in Parkin that result in a loss-of-function are the most frequent cause of familial Parkinson's disease and are responsible for a large percentage of autosomal recessive juvenile parkinsonism (AR-JP). In addition to its auto-ubiquitination, Parkin can also ubiquitinate various substrates, including the CDK2-interacting protein cyclin E1, the α -synuclein-interacting protein synphilin-1, the synaptic vesicle-associated CDC-rel1 and the p38 subunit of aminoacyl-tRNA synthetase/JTV-1. Ubiquitination of substrates occurs in conjunction with the E2

enzymes UbcH7, UbcH8 or UbcH13/Uev1.

Related Products Item # 23-021 UBE1, active, Item # 23-049 p38/JTV-1, Item # 23-047 UbcH7, active,

Item # 23-046 Parkin (c-Myc tagged), active

Selected References

Zhang Y. et al., Parkin Functions as an E2-Dependent Ubiquitin—Protein Ligase and Promotes the Degradation of the Synaptic Vesicle-Associated Protein, CDCrel-1. Proc Nat Acad Sci., 97: 13354–13359, 2000

Shimura H. *et al.*, Familial Parkinson Disease Gene Product, Parkin, is a Ubiquitin-Protein Ligase. Nature Genetics, *25*: 302-305, 2000

Imai Y. *et al.*, Parkin Suppresses Unfolded Protein Stress-Induced Cell Death through Its E3 Ubiquitin Protein Ligase Activity. J. Biol Chem., *275*: 35661-35664, 2000

Dev K.K. et al., Parkin-Associated Proteins and Parkinson's Disease. Neuropharmacology, 45: 1-13, 2003

Hampe C. *et al.*, Biochemical Analysis of Parkinson's Disease-Causing Variants of Parkin, an E3 Ubiquitin–Protein Ligase with Monoubiquitylation Capacity. Human Molecular Genetics, *15*: 2059-2075, 2006

Corti O. *et al.*, The p38 Subunit of the Aminoacyl-tRNA Synthetase Complex is a Parkin Substrate: Linking Protein Biosynthesis and Neurodegeneration. Hum Mol Genet., *12*: 1427–1437, 2003

Reviewed and approved by site quality representative.

Unless otherwise stated in our catalogue or other company documentation accompanying the product(s), our products are intended for research use only and are not to be used for any other purpose, which includes but is not limited to, unauthorized commercial uses, in vitro diagnostic uses, ex vivo or in vivo therapeutic uses or any type of consumption or application to humans or animals.

© 2014 Eurofins Pharma Discovery Services UK Limited is an independent member of Eurofins Discovery Services.