Camera & Optics

CX5 camera:

- · Scientific grade camera Made in Germany
- 5 megapixels resolution extendable to 20 megapixels
- · Passive cooling
- 16-bit 65 536 grey levels
- · USB-3 connection
- · Motorized zoom lens with feedback
- · Field of view: 26x21cm

Ease Of Use

One-Click-to-the-Image™ Auto-exposure Auto-focus Auto-lighting

Hardware

Smart Darkroom technology:

- · Software control of the lighting
- White light LED panels with automatic intensity adjustment
- UV cut-off filter and 7 positions filter wheel
- Slide-out PadBox with automatic recognition of the Application Pad inserted

Steel and stainless steel darkroom for long lasting robustness. Wide access door with UV safety shut-off

Software

Free Bio-Vision software for image acquisition with full GLP compliance. Molecular weight calculation, band quantification, colony counting, distance calculation, text annotation and image enhancement included.

CFR21 Part 11 ready

Applications

DNA and RNA gels and fluorescence stain imaging:

- Ethidium bromide, Sybr-Safe, Sybr-Green, Gel-Red, Gel-Green, Sybr-Gold, GFP, Pro-Q Emerald, Sypro ruby, FITC, DAPI
- Colorimetric stained protein gels, X-Ray film, autorads, SSCP gels, colony dish and flask imaging and other EPI white light applications
- Coomassie blue, Silver stain, Ponceau S Red, Copper stain...

Sky-Light Pad or Blue light conversion screen for DNA/RNA detection (to avoid damaging DNA):

· Sybr Safe, Sybr Green, eGFP

Technologies & Innovation

Apps Studio
3D Dynamic Scan
SuperResolution
ImageMaster™ assistant

Options

The Infinity is upgradeable to chemiluminescence

Several Application Pads can be easily inserted and removed inside the PadBox:

- UV Pad 312nm or dual wavelength 312 and 365nm
- Super-Bright Pad 312nm or dual wavelength 312 and 365nm
- Sky-Light Pad: LED blue light transilluminator, 470nm
- · White-Light Pad: white light LED transilluminator

White light or blue light conversion screen Advanced Bio-1D quantification software

CFR21 Part 11 module

