

TECHNICAL *Data Sheet*



environmental products for a greener planet

Parts Washer Regenerator Tabs

Key Benefits

- **Regenerates and Extends Solution Life:** Continuously digests oils, greases, and hydrocarbons to maintain bath cleanliness and extend service intervals
- **Improves System Efficiency:** Keeps filters, nozzles, and tanks cleaner for longer, ensuring steady flow and consistent cleaning power.
- **Reduces Odors and Organic Buildup:** Breaks down volatile organics and odor-causing residues for a cleaner, fresher work environment.
- **Supports Sustainable Operations:** Cuts chemical usage and fluid waste through natural biological cleaning and regeneration.
- **Safe & Easy to Use:** Non-toxic, non-pathogenic, non-corrosive microbes compatible with aqueous detergents and all metal surfaces.

Applications

- Industrial Parts Washers & Cleaning Systems
- Auto Repair & Maintenance Shops
- Machine & Manufacturing Facilities
- Fleet & Equipment Service Bays
- Oil Water Separators
- Fleet & Equipment Service Bays

Concentration (CFU/g):

- Tablets: 700 million CFU/gram
- Custom concentrations available upon request



Product Description

Parts Washer Regenerator Tabs are concentrated microbial tablets formulated to regenerate and extend the life of aqueous (water-based) parts washer solutions. This biological parts washer additive uses high-CFU, naturally occurring microbes that digest fats, oils, greases (FOGs), petroleum hydrocarbons, and other organic contaminants that accumulate during parts cleaning. By continuously breaking down these organics, the aqueous parts washer tablets keep wash fluid clean, reduce odors, and minimize costly fluid replacement and disposal.

Packaging Options

- Packs of 6 tablets
- Pails of 500 tablets
- Custom packaging available upon request

Shelf Life & Storage

- 2 years at room temperature
- 3 years refrigerated
- Store in a cool, dry location away from direct sunlight.
- Reseal container tightly after each use.

www.biogp.com

Application Guidelines

	Range	Optimal
pH Range	5.0 – 9.5	6.0 – 8.0
Temperature	45°F – 120°F (7°C – 49°C)	59°F – 104°F (15°C – 40°C)

Application Rates

Place 1 Parts Washer Regenerator Tab directly into the wash tank or fluid reservoir every 30 days for small to mid-size aqueous parts washers (10–40 gallons / 40–150 liters).

- For higher-load or continuous-use industrial parts washers, apply additional tablets or dose more frequently as needed to maintain bath clarity and consistent microbial activity.
- No mixing or agitation required, simply drop the tablet into the active wash bath.

Disclaimer: Results may vary depending on environmental conditions, application rates, and management practices. The manufacturer makes no guarantee of specific results. Seller's liability is limited to replacement of product or refund of purchase price. Manufacturer is not responsible for misuse, mishandling, or application under adverse conditions beyond its control. This product is not registered for pesticidal use with the U.S. Environmental Protection Agency. It is intended as a microbial inoculant only. Keep out of reach of children.

Contact Information

P: +1 (262) 285-4390

E: gogreen@biogp.com

520 Progress Drive
Belgium, WI 53004 USA