




INTRODUCING A SMARTER APPROACH TO ASSESSING

MALE INFERTILITY HEALTH

The most accurate, easy-to-use, and consistent
semen analysis result available today.




IMPROVE THE STANDARDIZATION OF YOUR SEMEN ANALYSIS PERFORMANCE

- Traditional Semen Analysis performance **varies from lab to lab¹**, even with optimal technique.
- **The best sample repeatability** due to superior field of view capability.
- **SEAMAN PRO™** reliably detects concentration samples versus other AI semen analysis systems.

System	Average CV (%)
SEAMAN PRO™	5.00
Other AI	14.03

DECREASED TURNAROUND TIME VERSUS TRADITIONAL SEMEN ANALYSIS

Expedite management and treatment earlier for fertility success, especially in IUI scenarios.



Results in
less than **5**
minutes

DEFINITIVE ANSWERS, CONFIDENT RESULTS

SEAMAN PRO™ provides definitive answers needed to ensure patients receive the **appropriate management and therapy as soon as possible**

Product	Description	Catalog #
Seaman PRO	System	#S-001
Seaman Slides	Fifty (50) Sperm Counting Slides	#CS-001

Reference:

1. Lemmens, Louise et al. "External quality control and training of semen analysis in the Netherlands: starting point for further reduction of outcome variability." Asian journal of andrology vol. 24,1 (2022): 15-20



6060 Center Dr. 10th Floor Suite 27
Playa Vista, CA 90045
www.checkcells.com

SEAMAN PRO™

CC-QA-2025-Rev 1