

Calf Serum

General information

Serum of the Newborn Calf category is obtained from 3 to 10 day old calves. In comparison, Calf Serum is obtained from animals which are not more than 12 months old.

Newborn Calf Serum and Calf Serum contain more immunoglobulins than Fetal Bovine Serum and have increased protein content. They are a cost effective alternative to Fetal Bovine Serum.

Newborn Calf Serum, Heat Inactivated has undergone heat inactivation, which involves heating the serum in a shaking water bath at exactly +56 °C for 30 minutes. This treatment inactivates the complement system, antibodies and other active enzymes.

Product specifications

pH	6.8 – 8.2
Osmolality	240 – 340 mOsm/kg
Endotoxin	As reported
Total Protein	5.0 – 8.5 g/dl
Albumin	As reported
Hemoglobin	As reported
Mycoplasma	Not detected
Virus tested for	PI-3, BVDV, BVDV-AB, BHV-1
Sterility	Tested
Storage	Store at ≤-15 °C.

Precautions and disclaimer

This product is for research use only. Please consult the Material Safety Data Sheet for information regarding hazards and safe handling practices.

Ordering information

Catalogue number	Pack size
LTC1000	1000 ml
LTCS500	500 ml
LTCS100	100 ml

Document version: v1.0. Copyright © 2026 LT Biotech, UAB. All rights reserved. Any redistribution or reproduction of part or all the content in any form is prohibited without the express written permission of LT Biotech, UAB.