



Revolute


Clarity · Speed · Confidence

Illuminate what matters

You bring the samples and ideas, we bring multiplex IF reagents to reveal the biological information you need
– spatial biology at the ease of standard IHC.

[revolute.de](https://www.revolute.de)

Workflow



What matters to you?



Clear Signal

Signal amplification enables bright staining for clear answers



Get results faster

Fast and simple protocol, staining within hours, not days



Maximal Ab flexibility

Use your own pre-validated Ab



Automatable

Reduces hands-on time; ensures reproducibility



Easy assay development

Develop assays within days, not months



Consistent reagent quality

Enables robust and consistent data generation

Don't worry about the workflow.

Focus on your science

Protect every tissue section

No antibody stripping.
No harsh treatments.

No conjugation

Use the antibodies
you already trust.
Continue to use
panels you already
established.

No rigid detection rules

Mix and match
primary antibodies.
Pair antibodies and
fluorophores however
you choose.

No lock-in

Run the protocol on
automated platforms
or manually. Image with
the instruments you
already have.

Get results faster

Parallel antibody staining and signal amplification allows you to perform experiments in half the time.

Staining time for 4-plex + nuclear stain [hrs]



Revolune



TSA

Available now

- 4-plex kit (secondary reagents for 4 biomarkers + nuclear counterstain)
- Suitable for rabbit antibodies

Available soon

- 5-, 6-, 7-plex kit
- Suitable for rabbit and mouse antibodies