

MRI Brain MS W & WO						
Technical Parameters						
Anatomic Coverage:	Skull base to vertex					
Coil:	Head					
Scanner Preference:	1.5T or 3T					
Contrast:	Pre-injection:			Post-injection:		
Scanning Plane:	Sagittal	Axial	Axial	Axial	Sagittal	Sagittal
Slice Thickness:	1.2mm	4mm	3mm	2mm	1mm	1mm
Gap:	0mm	1mm	0.3mm	-1mm	0mm	0mm
Weighting:	T1	T2	DWI	SWI	T2 FLAIR	T1
Mode:	3D TFE	2D TSE	2D EPI	3D GRE	3D IR-TSE	3D TSE
FAT SAT:	NO	NO	YES	NO	YES	YES
FOV:	24-25.6cm	23cm	23cm	23cm	23cm	23cm
MPR:	NO	NO	NO	NO	AXIAL, COR	AXIAL, COR
Notes:						