

MRI Brain Seizure W & WO							
Technical Parameters							
Anatomic Coverage:	Skull base to vertex						
Coil:	Head						
Scanner Preference:	1.5T or 3T						
Contrast:	Pre-injection:					Post-injection:	
	Scanning Plane:	Axial	Sagittal	Coronal Oblique	Coronal Oblique	Axial	Sagittal
Slice Thickness:	3mm	1.2mm	2mm	2.5mm	2mm	1mm	1mm
Gap:	0.3mm	0mm	0.2mm	0.25mm	-1mm	0mm	0mm
Weighting:	DWI	T1	T2 STIR	T2	SWI	T2 FLAIR	T1
Mode:	EPI	3D TFE	2D STIR TSE	2D TSE	3D GRE	3D IR-TSE	3D TSE
FAT SAT:	YES	NO	YES	NO	NO	YES	YES
FOV:	23cm	24-25.6cm	23cm	15cm	23cm	23cm	23cm
MPR:	NO	NO	NO	NO	NO	AXIAL, COR OBLIQ	AXIAL, COR OBLIQ
Notes:							