

	MRI Brain and IAC W & WO
--	--------------------------

Technical Parameters	
----------------------	--

Anatomic Coverage:	Skull base to vertex									
Coil:	Head									
Scanner Preference:	1.5 T or 3T									
Contrast:	Pre-Injection:							Post-Injection:		
	Scanning Plane:	Axial	Axial	Axial	Axial	Coronal	Sagittal	Axial (IACs)	Sagittal Brain	Axial Brain
Slice Thickness:	3mm	1mm	2mm	1mm	1mm	1mm	0.5mm	1mm	1mm	0.5mm
Gap:	0.3mm	0mm	-1mm	0mm	0mm	0mm	0mm	0mm	0mm	0mm
Weighting:	DWI	T2	SWI	T2	T2	T1	T2	T2 FLAIR	T1	T1
Mode:	EPI	3D TSE	3D GRE	STIR	STIR	3D TSE	3D TSE	3D IR-TSE	3D TFE	3D TSE
FAT SAT:	YES	NO	NO	YES	YES	NO	NO	YES	NO	YES
FOV:	23cm	23cm	23cm	23cm	23cm	23cm	17cm	23cm	23cm	17mm
MPR:	NO	NO	NO	NO	NO	AXIAL, COR	COR, SAG OBLIQ	AXIAL, COR	NO	COR, SAG
Notes:										