

MRI T-Spine WO			
Technical Parameters			
Anatomic Coverage:	C6-C7 to L1-L2		
Coil:	Torso and/or NV		
Scanner Preference:	1.5T or 3T		
Scanning Plane:	Sagittal	Sagittal	Sagittal
Slice Thickness:	1mm	3mm	3mm
Gap:	0mm	0.3mm	0.3mm
Weighting:	T2	T2	T1 FLAIR
Mode:	3D TSE	2D TSE	2D IR-TSE
FAT SAT:	NO	DIXON: IN PHASE	NO
FOV:			
MPR:	AXIAL 1/0	NO	NO
Notes:			