

MRI Ankle Right/Left WO

Technical Parameters

Anatomic Coverage:	Ankle					
Coil:	Ankle					
Scanner Preference:	1.5T or 3T					
Scanning Plane:	Axial	Axial	Sagittal	Sagittal	Coronal	Coronal
Slice Thickness:	3mm	3mm	3mm	3mm	3mm	3mm
Gap:	0.5mm	0.5mm	0.5mm	0.5mm	0.5mm	0.5mm
Weighting:	PD	PD FS	T1	STIR	PD	PD FS
Mode:	TSE	TSE	TSE	TSE	TSE	TSE
FAT SAT:	NO	YES	NO	YES	NO	YES
FOV:	15cm	15cm	15cm	15cm	15cm	15cm
MPR:	NO	NO	NO	NO	NO	NO
Notes:						

