

MRI Elbow Arthrogram			MRI Elbow WO (for reference)				
Technical Parameters			Technical Parameters				
Anatomic Coverage:	4cm above and below elbow joint		Anatomic Coverage:	4cm above and below elbow joint			
Coil:	Knee/Flex		Coil:	Knee/Flex			
Scanner Preference:	1.5T or 3T		Scanner Preference:	1.5T or 3T			
Contrast:	Post-injection:		Scanning Plane:	Axial	Axial	Coronal	Sagittal
			Scanning Plane:	Coronal	Sagittal	Slice Thickness:	3mm
Slice Thickness:	3mm	3mm	Gap:	0.5mm	0.5mm	0.5mm	0.5mm
Gap:	0.5mm	0.5mm	Weighting:	PD FS	T1	PD FS	PD FS
Weighting:	T1 FS	T1 FS	Mode:	TSE	TSE	TSE	TSE
Mode:	TSE	TSE	FAT SAT:	YES	NO	YES	YES
FAT SAT:	YES	YES	FOV:	14cm	14cm	14cm	14cm
FOV:	14cm	14cm	MPR:	NO	NO	NO	NO
MPR:	NO	NO	Notes:				
Notes:							

