

MRI Hip WO					
Technical Parameters					
Anatomic Coverage:	Above acetabulum to lesser trochanters				
Coil:	Torso				
Scanner Preference:	1.5T or 3T				
Scanning Plane:	Axial	Coronal	Coronal	Coronal	Sagittal
Slice Thickness:	4mm	5mm	4mm	4mm	4mm
Gap:	1mm	1mm	1mm	1mm	1mm
Weighting:	PD FS	STIR	T1	PD FS	PD FS
Mode:	TSE	TSE	TSE	TSE	TSE
FAT SAT:	YES	YES	NO	YES	YES
FOV:	18cm	Large FOV to include entire pelvis 30-40cm	18cm	18cm	18cm
MPR:	NO	NO	NO	NO	NO
Notes:					

