

MRI MSK Pelvis WO				
Technical Parameters				
Anatomic Coverage:	Iliac crest to lesser trochanters			
Coil:	Torso			
Scanner Preference:	1.5T or 3T			
Scanning Plane:	Axial	Coronal	Coronal	Sagittal
Slice Thickness:	4 mm	4 mm	4 mm	5 mm
Gap:	1mm	1mm	1mm	1mm
Weighting:	PDFS	T1	STIR	T2FS
Mode:	TSE	TSE	TSE	TSE
FAT SAT:	YES	NO	YES	YES
FOV:	36-40 cm	36-40 cm	36-40 cm	12-16 cm
MPR:	NO	NO	NO	NO
Notes:				

