

MRI Shoulder Arthrogram			MRI Shoulder WO (for reference)				
Technical Parameters			Technical Parameters				
Anatomic Coverage:	Entire shoulder including area of concern		Anatomic Coverage:	Entire shoulder			
Coil:	Shoulder		Coil:	Shoulder			
Scanner Preference:	1.5T or 3T		Scanner Preference:	1.5T or 3T			
Contrast:	Post-injection:						
			Scanning Plane:	Axial	Coronal	Sagittal	Sagittal
Scanning Plane:	Coronal	Sagittal	Slice Thickness:	3 mm	3 mm	3 mm	3 mm
Slice Thickness:	3 mm	3 mm	Gap:	0.5 mm	0.5 mm	0.5 mm	0.5 mm
Gap:	0.5 mm	0.5 mm	Weighting:	PD FS	PD FS	PD FS	T1
Weighting:	T1 FS	T1 FS					
			Mode:	TSE	TSE	TSE	TSE
Mode:	TSE	TSE					
			FAT SAT:	YES	YES	YES	NO
FAT SAT:	YES	YES					
			FOV:	16cm	16cm	16cm	16cm
FOV:	16cm	16cm					
MPR:	NO	NO					
MPR:	NO	NO	Notes:				
Notes:							