

MRI Shoulder WO				
Technical Parameters				
Anatomic Coverage:	Entire shoulder			
Coil:	Shoulder			
Scanner Preference:	1.5T or 3T			
Scanning Plane:	Axial	Coronal	Sagittal	Sagittal
Slice Thickness:	3 mm	3 mm	3 mm	3 mm
Gap:	0.5 mm	0.5 mm	0.5 mm	0.5 mm
Weighting:	PD FS	PD FS	PD FS	T1
Mode:	TSE	TSE	TSE	TSE
FAT SAT:	YES	YES	YES	NO
FOV:	16cm	16cm	16cm	16cm
MPR:	NO	NO	NO	NO
Notes:				