

MRI Tib/Fib WO				
Technical Parameters				
Anatomic Coverage:	Knees to Ankles			
Coil:	Torso			
Scanner Preference:	1.5T or 3T			
Scanning Plane:	Axial	Axial	Coronal	Sagittal
Slice Thickness:	5mm	5mm	5mm	5mm
Gap:	1mm	1mm	1mm	1mm
Weighting:	T1	PD FS	PD	STIR
Mode:	TSE	TSE	TSE	TSE
FAT SAT:	NO	YES	NO	YES
FOV:	40cm	40cm	40cm	40cm
MPR:	NO	NO	NO	NO
Notes:				

