

MRI Female Pelvis WO

Technical Parameters							
Anatomic Coverage:	Iliac crest to perinum						
Coil:	Torso						
Scanner Preference:	1.5T						
Scanning Plane:	Coronal	Axial	Axial	Sagittal	Coronal	Axial PRE	Axial
Slice Thickness:	5mm	5mm	5mm	5mm	5mm	5mm	5mm
Gap:	1mm	1mm	1mm	1mm	1mm	1mm	1mm
Weighting:	T2	T1	T2	T2	T2	T1	IN/OUT Phase
Mode:	SS-TSE BH	TSE	TSE	TSE	TSE FS	3D THRIVE	DWI/ADC
FAT SAT:	No	NO	No	NO	YES	YES	NO
FOV:	Top of kidneys to pelvis	20-24cm	20-24cm	24cm	20-24cm	Match AX TSE T2 Pelvis	Match AX TSE T2 Pelvis
MPR:							
Notes:							*If indication is for endometriosis or ovarian mass