

Minima *Pro*TM Stent System


Built to **Perform**. Built to **Reach Further**. Built for **Life**.

Redesigned ergonomic handle
with ease of use in mind

Now available in 10 mm balloon
same 4Fr profile across all 3 sizes

New hybrid open cell design
2.5x greater vessel retention*

70 cm working length
compatible with ≤ 65 cm long sheaths



Minima vs. Minima *Pro*™ Specifications

Specification	Minima	Minima <i>Pro</i> ™
Catheter effective length	63 cm	70 cm
Stent design	Closed-cell	Hybrid open-cell
Balloon options	6 mm, 8 mm	6 mm, 8 mm, 10 mm
Vessel contact points	~20	~50
Delivery handle	Twist Lock	Binary locking system
Packaging	Coiled	Straight
Profile	4 Fr	4 Fr
Initial implant range	5.1 – 8.5 mm	5.1 – 10.4 mm
Re-expansion to adult size	≤ 24 mm	≤ 24 mm
Guidewire Compatibility	.014" and .018"	.014" and .018"
FDA approved indication	CoA & BPAS	CoA & BPAS

Ordering Information

Part Number	FG-0009	FG-0010	FG-0011
Stent Diameter x Length	6 mm x 16 mm	8 mm x 16 mm	10 mm x 16 mm
Nominal Pressure	8 atm	10 atm	8 atm
Rated Burst Pressure	14 atm	14 atm	10 atm

Important Information:

Indications, contraindications, warnings, and instructions for use can be found with the product labeling and packaging materials supplied with each device or at renatamedical.com/manuals. Federal law restricts this product for sale by or on the order of a physician. This product is intended only for use by or under the direction of a physician and approved for U.S. use only.

Illustrations are artistic renderings only and should not be considered as engineering drawings or photographs. Photos are on file at Renata Medical, Inc., 4675 MacArthur Court, Suite 1150, Newport Beach, CA 92660.

Minima, Renata Medical, and its logo are registered trademarks of Renata Medical, Inc.

*Data on File (TS-0018 Rev. 0).

MKT-0018 Rev. 0

