MYCO MEDICAL
TheSharpChoice
NMSDC Certified Diversity Supplier

Reli®

EPIDURAL Needles

RELI® Epidural Needles are the value alternative to leading brands and offer superior quality in a broad range of sizes.




Fixed wing family (TUFW)

- Needle gauges available: 16G to 22G
- Needle lengths available: 2.5" to 6"
- Precision point-grinding for comfort and non-coring benefit
- Wide metric marks for maximum visibility and accuracy
- Available with plastic or metal stylet with integrated wings or detachable wings
- Large, clear hub for enhanced tactile feel
- Color-coded stylet caps
- Available in sterile and non-sterile packaging


16G White	17G Violet	18G Pink	20G Yellow	22G Black

BRAND (COMPANY)	CLEAR HUBS	COLOR-CODED STYLETS	FITTED, CLOSE TOLERANCE STYLET	INTEGRATED & DETACHABLE WING OPTIONS
Reli® Epidural Needles (MYCO Medical)	✓	✓	✓	✓
Perisafe and Weiss (BD)		✓	✓	✓
Perifix (B. Braun)	✓	✓	✓	
Tuohy Epidural (Avanos)	✓	✓	✓	

1.919.460.2535 www.mycomedical.com



TU18G detachable wing style, shown without wing



Precision point-grinding



TU detachable wing style, shown with wings

ORDER INFORMATION

Unit of measure · **Sterile:** 25 ea/bx; 4 bx/cs · **Non-Sterile:** 100 ea/bx

WITH DETACHABLE WING

GAUGE	LENGTH	STERILE MYCO Code	NON-STERILE MYCO Code
16G	3.5"	TU16G351	TU16G350
17G	3.5"	TU17G351	TU17G350
17G	6"	TU17G601	TU17G600
18G	2.5"	TU18G251	TU18G250
18G	3.5"	TU18G351	TU18G350
18G	4.5"	TU18G451	TU18G450
18G	6"	TU18G601	TU18G600
20G	3.5"	TU20G351	TU20G350
20G	4.5"	TU20G451	TU20G450
20G	6"	TU20G601	TU20G600
22G	2.5"	TU22G251	TU22G250
22G	3.5"	TU22G351	TU22G350

WITH INTEGRATED WING

GAUGE	LENGTH	STERILE MYCO Code	NON-STERILE MYCO Code
17G	3.5"	TUFW17G351	TUFW17G350
18G	3.5"	TUFW18G351	TUFW18G350
20G	3.5"	TUFW20G351	TUFW20G350

MYCO MEDICAL
TheSharpChoice
NMSDC Certified Diversity Supplier

1.919.460.2535
2015 Production Drive
Apex, NC 27539 USA
www.mycomedical.com
info@mycomedical.com

Scan to learn more
or receive complimentary evaluation samples
https://mycomed.co/epidural-needles

Some medical devices are not licensed for sale in Canada. © 2026 MYCO Medical, Inc. All rights reserved.
MYCO® is a registered trademark of MYCO Medical, Inc. · LIT-00195 Rev 3 Eff 11.04.2026