VIVIT SUPPLY

MEDICAL FACE MASK EARLOOPS • BLUE ASTM LEVEL 1



101119EM

Combining comfort and breathability, this ASTM Level 1 mask offers standard fluid resistance and filtration efficiency to keep clinicians protected during procedures and patient care. This mask can also be worn when experiencing respiratory symptoms, such as coughing or sneezing, to help limit the spread of droplets containing pathogens.

BFE	PFE	DELTA P $\Delta P(mmH_2O/cm^2)$	FLUID RESISTANCE	FLAME SPREAD
> 95%	> 80%	< 5.0	80 mm Hg	Class 1

MASKS/BOX	BOXES/CASE	MASKS/CASE
50	10	500

① Ultra-soft Earloop

- Earloop material is soft and smooth for maximum comfort.
- Round design makes this mask one of the most comfortable options available.

② High Quality Nosepiece

- Aluminum nosepiece forms smoothly across the nose and cheekbones.
- Remains securely in place and reduces fogging of eyewear.

③ Responsive Inner Layer Material

- Padded feel for ultimate comfort.
- Flexible material that moves with you.

NON STERILE · 2 · LATEX Not made with natural rubber latex

MADE IN CANADA

VIVIT Supply LLC
200 Continental Dr., Suite 401
Newark, DE 19713 USA

W vivitsupply.com

*Registered Trademark or Trademark of PRIMED Medical Products Inc. or its affiliates.