VIVIT* SUPPLY

MEDICAL FACE MASK EARLOOPS • OPAL ASTM LEVEL 3



301149EM

Using only premium quality materials, this ASTM Level 3 mask offers the highest levels of fluid resistance and filtration efficiency. Rigorously tested against the ASTM F2100 Standards, clinicians can be confident they are protected during medical procedures and patient care in areas such as the ER and ICU.

BFE	PFE	DELTA P $\Delta P(mmH_2O/cm^2)$	FLUID RESISTANCE	FLAME SPREAD
> 98%	> 85%	< 6.0	160 mm Hg	Class 1

MASKS/BOX	BOXES/CASE	MASKS/CASE
50	10	500

① Ultra-soft Earloop

- Earloop material is soft and smooth for maximum comfort.
- Round design makes this mask one of the most comfortable options available.

② High Quality Nosepiece

- Aluminum nosepiece forms smoothly across the nose and cheekbones.
- Remains securely in place and reduces fogging of eyewear.

③ Smooth Naturally Breathable Inner Layer

- Provides a light and cool feel to the mask, reducing chances of irritation throughout extended wear.
- Inner layer maintains its structure, allowing the mouth to feel open and airy.

NON STERILE · 2 · LATEX Not made with natural rubber latex

MADE IN CANADA

VIVIT Supply LLC
200 Continental Dr., Suite 401
Newark, DE 19713 USA

W vivitsupply.com

*Registered Trademark or Trademark of PRIMED Medical Products Inc. or its affiliates.