VIVIT SUPPLY

MEDICAL FACE MASK
EARLOOPS • SOFT • BLUE
ASTM LEVEL 3



30211EM

As a physical barrier against contaminants, this ASTM Level 3 mask offers the highest level of fluid resistance and filtration efficiency to keep clinicians protected during surgical procedures. The premium inner layer—a soft, irritant-free material—makes this the perfect mask for individuals with skin sensitivities.

BFE	PFE	DELTA P $\Delta P(mmH_2O/cm^2)$	FLUID RESISTANCE	FLAME SPREAD
> 98%	> 85%	< 6.0	160 mm Hg	Class 1

MASKS/BOX	BOXES/CASE	MASKS/CASE
50	10	500

① Ultra-soft Earloop

- Earloop material is soft and smooth for maximum comfort.
- Round design makes this mask one of the most comfortable options available.

② High Quality Nosepiece

- Aluminum nosepiece forms smoothly across the nose and cheekbones.
- Remains securely in place and reduces fogging of eyewear.

③ Soft Irritant Free Inner Layer

- Inner layer has a silky finish that is delicate yet comfortable, making it ideal for those with skin sensitivities.

NON STERILE · 2 · LATEX Not made with natural rubber latex

VIVIT Supply LLC
200 Continental Dr., Suite 401
Newark, DE 19713 USA

W vivitsupply.com

*Registered Trademark or Trademark of PRIMED Medical Products Inc. or its affiliates.