VIVIT* SUPPLY

MEDICAL FACE MASK
EARLOOPS • ANTI-FOG FOAM • OPAL
ASTM LEVEL 3



313149EM

This ASTM Level 3 mask offers the highest levels of fluid resistance and filtration efficiency with the added benefit of an anti-fog foam strip for enhanced visibility when worn with eyewear. Clinicians reach for this mask when performing medical procedures with a high risk of fluid splash in areas including the ER, ICU, and trauma.

BFE	PFE	DELTA P $\Delta P(mmH_2O/cm^2)$	FLUID RESISTANCE	FLAME SPREAD
> 98%	> 85%	< 6.0	160 mm Hg	Class 1

MASKS/BOX	BOXES/CASE	MASKS/CASE
50	10	500

1 Ultra-soft Earloop

- Earloop material is soft and smooth for maximum comfort.
- Round design makes this mask one of the most comfortable options available.

2 High Quality Nosepiece with Ultra Soft Foam

- Aluminum nosepiece and soft foam strip work together to form smoothly across the nose and cheekbones.
- Remains securely in place to create a custom and secure seal providing enhanced resistance to fogging of eyewear and a cushion-like feel to the mask.

3 Smooth Naturally Breathable Inner Layer

- Provides a light and cool feel to the mask, reducing chances of irritation throughout extended wear.
- Inner layer maintains its structure, allowing the mouth to feel open and airy.

NON STERILE · 2 · LATEX Not made with natural rubber latex

MADE IN CANADA

VIVIT Supply LLC
200 Continental Dr., Suite 401
Newark, DE 19713 USA

W vivitsupply.com

*Registered Trademark or Trademark of PRIMED Medical Products Inc. or its affiliates.