VIVIT SUPPLY

MEDICAL FACE MASK
TIES • DUCKBILL • ANTI-FOG FOAM • OPAL
ASTM LEVEL 3



31334TM

Rigorously tested against the ASTM F2100 Standards, this ASTM Level 3 mask offers maximum protection for both fluid resistance and filtration efficiency. Clinicians can be confident when using this mask during surgical procedures including trauma, ER, and orthopedic cases. Lightweight materials, combined with a unique duckbill shape, create a large breathing chamber allowing clinicians to experience maximum breathability.

BFE	PFE	DELTA P $\Delta P(mmH_2O/cm^2)$	FLUID RESISTANCE	FLAME SPREAD
> 98%	> 85%	< 6.0	160 mm Hg	Class 1

MASKS/BOX	BOXES/CASE	MASKS/CASE
50	6	300

1 Premium Tie Material and Design

- Tie material that is comfortable on the skin throughout extended wear.
- Long ties create a flexible fit while tying over top of a bouffant or surgical cap.

2 High Quality Nosepiece with Ultra Soft Foam

- Aluminum nosepiece and soft foam strip work together to form smoothly across the nose and cheekbones.
- Remains securely in place to create a custom and secure seal providing enhanced resistance to fogging of eyewear and a cushion-like feel to the mask.

3 Optimized Duckbill Chamber

- Large breathing chamber to keep the mask away from the nose and mouth during use.
- Provides maximum comfort, and breathability ideal for surgical procedures and extended wear.

NON STERILE · LATEX Not made with natural rubber latex

VIVIT Supply LLC
200 Continental Dr., Suite 401
Newark, DE 19713 USA

W vivitsupply.com

*Registered Trademark or Trademark of PRIMED Medical Products Inc. or its affiliates.