VIVIT SUPPLY

MEDICAL FACE MASK
TIES • ANTI-FOG TAPE • GREEN
ASTM LEVEL 3



33313TM

For use in environments with high fluid exposure, this ASTM Level 3 mask establishes a physical barrier between clinicians and contaminants from their surrounding environment. The pressure-sensitive adhesive strip allows for a tailored, secure seal to prevent the fogging of eyewear, providing a clear line of sight during surgical procedures.

BFE	PFE	DELTA P $\Delta P(mmH_2O/cm^2)$	FLUID RESISTANCE	FLAME SPREAD
> 98%	> 85%	< 6.0	160 mm Hg	Class 1

MASKS/BOX	BOXES/CASE	MASKS/CASE
50	6	300

① Premium Tie Material and Design

- Tie material that is comfortable on the skin throughout extended wear.
- Long ties create a flexible fit while tying over top of a bouffant or surgical cap.

② High Quality Nosepiece with Anti-fog Tape

- Nosepiece combined with a skin-friendly adhesive strip creates the ultimate anti-fog face mask, ensuring a complete seal across the nose and cheekbones to prevent the fogging of eyewear.

③ Smooth Naturally Breathable Inner Layer

- Provides a light and cool feel to the mask, reducing chances of irritation throughout extended wear.
- Inner layer maintains its structure, allowing the mouth to feel open and airy.

NON STERILE · LATEX Not made with natural rubber latex

VIVIT Supply LLC
200 Continental Dr., Suite 401
Newark, DE 19713 USA

W vivitsupply.com

*Registered Trademark or Trademark of PRIMED Medical Products Inc. or its affiliates.