VIVIT SUPPLY

EYEWEAR
FRAME • LENS • ASSEMBLED
MULTICOLORED



11EA / 11ED/ 11EF / 11EL

VIVIT Supply eyewear is a convenient and cost-effective solution for clinicians, protecting vulnerable mucous membranes from exposure to harmful fluids and particles that spread infection and illness. With pre-assembled options, reusable frames and single-use lenses, the VIVIT Supply eyewear line is well-equipped for a variety of procedures with a risk of fluid splash. The reusable wrap-around frames provide a secure fit that is both comfortable and accommodating when worn with prescription glasses, and the lenses come with double-sided anti-fog coating help maintain an unobstructed view during medical procedures.

SKU#	PRODUCT TYPE	PIECES/BOX	BAGS/CASE	PIECES/CASE
11EA	Assembled	10	5	50
11ED	Dispenser	25	4	100
11EF	Frame	10	10	100
11EL	Lenses	25	10	250

① Premium Anti-Fog Lenses

- Increased optical clarity over untreated lenses for a clear line of sight, unencumbered by moisture buildup.
- Protection from direct and peripheral spray of fluids with wide lens coverage.

② Comfortable and Secure Fit

- Comfortable and secure fit, large enough to wear over perscription glasses.

③ Colored Frames

- Easily identifiable, reusable frame that reduces wastage.
- Multicolored:

NON STERILE · LATEX Not made with natural rubber latex

VIVIT Supply LLC
200 Continental Dr., Suite 401
Newark, DE 19713 USA

W vivitsupply.com