
## Setting up Dexcom Clarity


# How can I share my sensor data with my healthcare team?


You can now easily share your data with your healthcare team through the Dexcom sensor app or a web browser. Obtain the Clarity Clinic Code from your healthcare provider and follow these simple steps.


### Within the Dexcom G7 App:

1. Open your Dexcom sensor app on your smartphone
2. Tap the Connections Tab
3. Tap Clarity Clinic and tap “Continue”
4. Enter the Clarity Clinic Code when prompted.
5. Clinic ID: 5716343636
6. Confirm clinic details and tap “Done” when finished.
7. Clarity clinic will remain an active connection in the Connections Tab of your Dexcom CGM App and you can use it to stop sharing and/or share with other Clarity clinics.





9:41
Connections
Sensor
4 days left
Clarity Clinic
Share data using Dexcom Clarity.
Available connections
Apple Health
Clarity Clinic
Clarity Clinic lets you share your CGM data with your clinic.
Tap Continue to start sharing with your clinic and to see your connected clinics.
Continue



9:41
Done
Clarity Clinic
Enter clinic code
Enter the clinic code to begin sharing CGM data.
Clinic code
Continue
Ask your clinic for their clinic code.
Connected clinics
Caring Touch Center




*Sample screens from Dexcom G7 app*

You can stop sharing your data anytime by selecting the clinic under Connected Clinics and tap Stop Sharing.






If your health care provider does not have a Dexcom Clarity clinic account you can download or print PDF reports through Dexcom Clarity.

## Through a web browser

1. Navigate to connect.dexcom.com (.jp for Japan and Singapore; .eu for ROW)
2. Login with your Dexcom Credentials.
3. Enter the Clarity Clinic Code when prompted.
4. Confirm clinic details and tap “Done” when finished.
5. Clarity clinic will remain an active connection in the Connections Tab of your Dexcom CGM App and you can use it to stop sharing and/or share with other Clarity clinics.

**This applies to both Dexcom G7 and Dexcom G7 15 Day CGM Systems.**